



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### Blood DNA Purification Kit (Pre-filled)

**Catalog Number:** U0071

**Regulatory Status:** For research use only (RUO)

**Product Description:** Blood DNA Purification Kit (Pre-filled) is designed for the optimal automated extraction of DNA from fresh human whole blood and buffy coat.

**Applications:** Purify  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Instruction:** Components under room temperature (15~35 °C) can be stored until the expiration date labeled on the box. The Proteinase K was transported at room temperature. When received, please store at -20 °C. Repeating of freezing and thawing may cause the activity decay of Proteinase K.