



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### Heat Pretreatment Solution (EDTA)

**Catalog Number:** U0542

**Regulation Status:** For research use only (RUO)

**Product Description:** The Heat Pretreatment Solution (EDTA) is intended to be used to heat pretreatment of formalin-fixed, paraffin-embedded specimens in CISH.

**Regulatory Status:** For research use only (RUO)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** Chromogenic *In Situ* Hybridization (FFPE Tissue)

**Storage Instruction:** Store at 2-8°C in an upright position. Return to storage conditions immediately after use.