



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### DAPI solution

**Catalog Number:** U0553

**Regulation Status:** For research use only (RUO)

**Product Description:** DAPI solution is used for dead cell staining. DAPI is utilized for the detection of mitochondrial DNA in yeast, chloroplast DNA, virus DNA, micoplasm DNA and chromosomal DNA. The excitation and emission wavelengths of DAPI-DNA complex are 360 nm and 460 nm, respectively.

**Theoretical MW (kDa):** 350.25

**Regulation Status:** For research use only (RUO)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Storage Instruction:** Store at 0-5°C, Protect from light.