



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

Hoechst 33342 solution

Catalog Number: U0556

Regulation Status: For research use only (RUO)

Product Description: Hoechst 33342 solution is used as nuclear dye. Hoechst dye is cell membrane permeable and stain DNA to emit intense blue fluorescence. The excitation and emission wavelengths of Hoechst-DNA complex are 350 nm and 460 nm, respectively.

Theoretical MW (kDa): 561.9299999999999

Regulation Status: For research use only (RUO)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Storage Instruction: Store at 0-5°C, Protect from light.