



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

Human CD3/CD28 ActiveBeads™

Catalog Number: U0576

Regulatory Status: For research use only (RUO)

Product Description: Human CD3/CD28 ActiveBeads™
induce T cell activation in enriched T cells or PBMCs.

Suitable Sample: Enriched T cells or PBMCs

Cell Type: T Cells

Product Type: CD3/CD28 T-Activator

Particle Core: Magnetic beads

Antibody Conjugation: Humanized anti-human CD3
and anti-human CD28 monoclonal antibodies

Capacity: Up to 1×10^8 / enriched T cells or up to 2×10^8
PBMCs

Regulatory Status: For research use only (RUO)

Concentration: 4×10^8 Beads/ mL

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Storage Buffer: In PBS containing Human Serum
Albumin (HSA), pH 7.4.

Storage Instruction: Store at 4°C. Do not freeze.