



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

GST agarose beads

Catalog Number: U0580

Regulation Status: For research use only (RUO)

Product Description: GST agarose beads is based on anti-GST covalently coupled to the surface of agarose beads.

Ligand: Anti-GST single domain antibody fragment (VHH, Nanobody)

Bead Size: ~40 um

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: The GST agarose beads recognizes GST selectively.

Form: Liquid

Recommend Usage: Co-Immunoprecipitation

Enzyme activity

Immunoprecipitation

Mass Cytometry

The optimal working dilution should be determined by the end user.

Storage Buffer: 50% slurry in PBS (20% ethanol).

Storage Instruction: Store at 4°C. Do not freeze.