



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### TurboGFP agarose beads

**Catalog Number:** U0584

**Regulation Status:** For research use only (RUO)

**Product Description:** TurboGFP agarose beads is based on anti-TurboGFP covalently coupled to the surface of agarose beads.

**Ligand:** Anti-TurboGFP single domain antibody fragment (VHH, Nanobody)

**Bead Size:** ~40 um

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** The TurboGFP agarose beads recognizes TurboGFP selectively.

**Form:** Liquid

**Recommend Usage:** Co-Immunoprecipitation  
Enzyme activity  
Immunoprecipitation  
Mass Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** 50% slurry in PBS (20% ethanol).

**Storage Instruction:** Store at 4°C. Do not freeze.