



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

BB19012 : Prostate intraepithelial neoplasia and cancer tissue array, including TNM, clinical stage and pathology grade, 24 cases/48 cores

Microarray Panel	Prostate intraepithelial neoplasia and carcinoma tissue microarray, containing 21 cases of prostate intraepithelial neoplasia and 3 prostate adenocarcinoma, duplicate cores per case		
Cores	48	BB19012 H&E 	
Cases	24		
Row number	6		
Column number	8		
Core Diameter (mm)	2		
Thickness (µm)	5		
Tissue Array Type	FFPE		
Species	Human		
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.		
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, “*” means invalid core; “-” means no applicable or negative in IHC markers.		

	1	2	3	4	5	6	7	8
A	Pro	Pro	Pro	Pro	Pro	Pro	Pro	Pro
B	Pro	Pro	Pro	Pro	Pro	Pro	Pro	Pro
C	Pro	Pro	Pro	Pro	Pro	Pro	Pro	Pro
D	Pro	Pro	Pro	Pro	Pro	Pro	Pro	Pro
E	Pro	Pro	Pro	Pro	Pro	Pro	Pro	Pro
F	Pro	Pro	Pro	Pro	Pro	Pro	Pro	Adr

Legend: Pro - Prostate

● - Malignant tumor, ● - Malignant tumor (stage II), ● - Malignant tumor (stage III), ● - PIN

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	70	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
A2	2	70	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
A3	3	76	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
A4	4	76	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
A5	5	68	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
A6	6	68	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
A7	7	80	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
A8	8	80	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
B1	9	70	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
B2	10	70	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
B3	11	77	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
B4	12	77	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
B5	13	71	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
B6	14	71	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
B7	15	69	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
B8	16	69	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
C1	17	67	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
C2	18	67	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	

C3	19	70	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
C4	20	70	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
C5	21	63	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
C6	22	63	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
C7	23	79	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
C8	24	79	M	Prostate	Prostate intraepithelial neoplasia	-	1	-	PIN	
D1	25	70	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
D2	26	70	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
D3	27	55	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
D4	28	55	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
D5	29	74	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
D6	30	74	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
D7	31	78	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
D8	32	78	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
E1	33	70	M	Prostate	Prostate intraepithelial neoplasia	-	3	-	PIN	
E2	34	70	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
E3	35	71	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
E4	36	71	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
E5	37	81	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
E6	38	81	M	Prostate	Prostate intraepithelial neoplasia	-	2	-	PIN	
E7	39	75	M	Prostate	Prostate intraepithelial neoplasia	-	3	-	PIN	
E8	40	75	M	Prostate	Prostate intraepithelial neoplasia	-	3	-	PIN	
F1	41	64	M	Prostate	Prostate intraepithelial neoplasia	-	3	-	PIN	
F2	42	64	M	Prostate	Prostate intraepithelial neoplasia	-	3	-	PIN	
F3	43	71	M	Prostate	Adenocarcinoma	T2N0M0	1	II	malignant	
F4	44	71	M	Prostate	Adenocarcinoma	T2N0M0	1	II	malignant	
F5	45	73	M	Prostate	Adenocarcinoma	T3N0M0	1	III	malignant	
F6	46	73	M	Prostate	Adenocarcinoma	T3N0M0	1	III	malignant	
F7	47	82	M	Prostate	Adenocarcinoma	T2N0M0	2	II	malignant	
F8	48	82	M	Prostate	Adenocarcinoma	T2N0M0	1	II	malignant	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada

1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU

Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World

184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland

Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland

T: +49 (0) 69 779099
F: +49 (0) 69 13376880