



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

LV243 : Liver cancer with liver tissue array, including pathology grade, TNM and clinical stage, 12 cases/24 cores

Microarray Panel	Hepatocellular carcinoma with liver tissue microarray, containing 10 cases of hepatocellular carcinoma and 2 liver tissues, duplicated cores per case		
Cores	24	LV243 H&E 	
Cases	12		
Row number	4		
Column number	6		
Core Diameter (mm)	1.5		
Thickness (µm)	5		
QA/QC	H&E, IHC		
Tissue Array Type	FFPE		
Species	Human		
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.		
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheets, “*” means invalid core; “-” means no applicable or negative in IHC markers.		

	1	2	3	4	5	6
A	Liv	Liv	Liv	Liv	Liv	Liv
B	Liv	Liv	Liv	Liv	Liv	Liv
C	Liv	Liv	Liv	Liv	Liv	Liv
D	Liv	Liv	Liv	Liv	Liv	Sk1

Legend: Liv - Liver

● - Malignant tumor,
 ● - Malignant tumor (stage II),
 ● - Malignant tumor (stage IIIA),
 ● - Malignant tumor (stage IVA),
 ● - Normal tissue

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Notes	Image
A1	1	47	M	Liver	Hepatocellular carcinoma	T3N0M0	2	IIIA	Malignant		
A2	2	47	M	Liver	Hepatocellular carcinoma	T3N0M0	2	IIIA	Malignant		
A3	3	49	M	Liver	Hepatocellular carcinoma	T2N0M0	2	II	Malignant		
A4	4	49	M	Liver	Hepatocellular carcinoma	T2N0M0	2	II	Malignant		
A5	5	31	M	Liver	Hepatocellular carcinoma	T3N0M0	2	IIIA	Malignant		
A6	6	31	M	Liver	Hepatocellular carcinoma	T3N0M0	2	IIIA	Malignant		
B1	7	62	M	Liver	Hepatocellular carcinoma	T3N0M0	2	IIIA	Malignant		
B2	8	62	M	Liver	Hepatocellular carcinoma	T3N0M0	2	IIIA	Malignant		
B3	9	38	M	Liver	Hepatocellular carcinoma	T3N0M0	2	IIIA	Malignant		
B4	10	38	M	Liver	Hepatocellular carcinoma	T3N0M0	2	IIIA	Malignant		
B5	11	61	M	Liver	Hepatocellular carcinoma	T2N0M0	2	II	Malignant		
B6	12	61	M	Liver	Hepatocellular carcinoma	T2N0M0	2	II	Malignant		
C1	13	59	M	Liver	Hepatocellular carcinoma	T2N0M0	2	II	Malignant		
C2	14	59	M	Liver	Hepatocellular carcinoma	T2N0M0	2	II	Malignant		
C3	15	43	M	Liver	Hepatocellular carcinoma	T2N0M0	3	II	Malignant		
C4	16	43	M	Liver	Hepatocellular carcinoma	T2N0M0	3	II	Malignant		
C5	17	48	M	Liver	Hepatocellular carcinoma	T3N1M0	3	IVA	Malignant		
C6	18	48	M	Liver	Hepatocellular carcinoma	T3N1M0	3	IVA	Malignant		
D1	19	78	F	Liver	Hepatocellular carcinoma	T3N0M0	3	IIIA	Malignant		
D2	20	78	F	Liver	Hepatocellular carcinoma	T3N0M0	3	IIIA	Malignant		
D3	21	16	M	Liver	Liver tissue	-	-	-	Normal		
D4	22	16	M	Liver	Liver tissue	-	-	-	Normal		
D5	23	40	M	Liver	Liver tissue	-	-	-	Normal		
D6	24	40	M	Liver	Liver tissue	-	-	-	Normal		
-	0	58	M	Skin	Malignant melanoma (tissue marker)	-	-	-	Malignant	-	