



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

LV486 : Liver cancer tissue array with matched adjacent tumor tissue, including TNM and pathology grade, 16 cases/48 cores

Microarray Panel	Liver hepatocellular carcinoma and matched adjacent tissue microarray, containing 16 cases of liver carcinoma and matched adjacent tissue, triplicate cores per case	
Cores	48	LV486 H&E 
Cases	16	
Row number	6	
Column number	8	
Core Diameter (mm)	2	
Thickness (µm)	5	
Quality Control	Anti-Actin confirmed	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, “**” means invalid core; “-” means no applicable or negative in IHC markers.	

	1	2	3	4	5	6	7	8
A	Liv	Liv	Liv	Liv	Liv	Liv	Liv	Liv
B	Liv	Liv	Liv	Liv	Liv	Liv	Liv	Liv
C	Liv	Liv	Liv	Liv	Liv	Liv	Liv	Liv
D	Liv	Liv	Liv	Liv	Liv	Liv	Liv	Liv
E	Liv	Liv	Liv	Liv	Liv	Liv	Liv	Liv
F	Liv	Liv	Liv	Liv	Liv	Liv	Liv	Liv

Legend: Liv - Liver

● - AT, ● - Malignant tumor, ● - Malignant tumor (stage II), ● - Malignant tumor (stage III), ● - Malignant tumor (stage IIIc)

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	46	M	Liver	Hepatocellular carcinoma	T2N0M0	1	II	malignant	
A2	2	49	M	Liver	Hepatocellular carcinoma	T2N0M0	2	II	malignant	
A3	3	50	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
A4	4	71	M	Liver	Hepatocellular carcinoma	T3N0M0	1	III	malignant	
A5	5	69	F	Liver	Hepatocellular carcinoma with necrosis	T4N1M0	2	IIIc	malignant	
A6	6	37	F	Liver	Hepatocellular carcinoma	T3N0M0	1	III	malignant	
A7	7	56	M	Liver	Hepatocellular carcinoma	T3N0M0	3	III	malignant	
A8	8	60	M	Liver	Hepatocellular carcinoma	T2N0M0	1	II	malignant	
B1	9	46	M	Liver	Hepatocellular carcinoma	T2N0M0	1	II	malignant	
B2	10	49	M	Liver	Hepatocellular carcinoma	T2N0M0	2	II	malignant	
B3	11	50	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
B4	12	71	M	Liver	Hepatocellular carcinoma	T3N0M0	1	III	malignant	
B5	13	69	F	Liver	Hepatocellular carcinoma with necrosis (sp... e)	T4N1M0	2	IIIc	malignant	

B6	14	37	F	Liver	Hepatocellular carcinoma	T3N0M0	1	III	malignant	
B7	15	56	M	Liver	Hepatocellular carcinoma	T3N0M0	3	III	malignant	
B8	16	60	M	Liver	Hepatocellular carcinoma	T2N0M0	1	II	malignant	
C1	17	46	M	Liver	Cancer adjacent tissue (chronic hepatitis)	-	-	-	AT	
C2	18	49	M	Liver	Cancer adjacent tissue (liver cirrhosis)	-	-	-	AT	
C3	19	50	M	Liver	Cancer adjacent tissue (hepatic tissue)	-	-	-	AT	
C4	20	71	M	Liver	Cancer adjacent tissue (hepatic tissue)	-	-	-	AT	
C5	21	69	F	Liver	Cancer adjacent tissue (hepatic tissue)	-	-	-	AT	
C6	22	37	F	Liver	Cancer adjacent tissue (hepatic tissue)	-	-	-	AT	
C7	23	56	M	Liver	Cancer adjacent tissue (liver cirrhosis sparse)	-	-	-	AT	
C8	24	60	M	Liver	Cancer adjacent tissue (liver cirrhosis)	-	-	-	AT	
D1	25	63	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
D2	26	43	F	Liver	Hepatocellular carcinoma	T3N0M0	3	III	malignant	
D3	27	56	M	Liver	Hepatocellular carcinoma	T2N0M0	1	II	malignant	
D4	28	52	F	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
D5	29	67	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
D6	30	64	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
D7	31	35	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
D8	32	43	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
E1	33	63	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
E2	34	43	F	Liver	Hepatocellular carcinoma	T3N0M0	3	III	malignant	
E3	35	56	M	Liver	Hepatocellular carcinoma	T2N0M0	1	II	malignant	
E4	36	52	F	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
E5	37	67	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
E6	38	64	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
E7	39	35	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
E8	40	43	M	Liver	Hepatocellular carcinoma	T3N0M0	2	III	malignant	
F1	41	63	M	Liver	Cancer adjacent tissue (hepatic tissue sparse)	-	-	-	AT	
F2	42	43	F	Liver	Cancer adjacent tissue (liver cirrhosis with fatty degeneration)	-	-	-	AT	
F3	43	56	M	Liver	Cancer adjacent tissue (chronic hepatitis)	-	-	-	AT	
F4	44	52	F	Liver	Cancer adjacent tissue (hepatic tissue)	-	-	-	AT	
F5	45	67	M	Liver	Cancer adjacent tissue (liver cirrhosis)	-	-	-	AT	
F6	46	64	M	Liver	Cancer adjacent tissue (liver cirrhosis)	-	-	-	AT	
F7	47	35	M	Liver	Cancer adjacent tissue (chronic hepatitis)	-	-	-	AT	
F8	48	43	M	Liver	Cancer adjacent tissue (liver cirrhosis)	-	-	-	AT	
-	0	58	M	Skin	Malignant melanoma (tissue marker)		-		Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
 184 Park Drive, Milton Park
 Abingdon OX14 4SE.
 T: +44 (0) 1235 828 200
 F: +44 (0) 1235 820 482

North America 
 1035 Cambridge Street,
 Cambridge, MA 02141.
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

Europe 
 Berenkoog 41,
 1822 BH Alkmaar,
 Netherlands
 T: +31 (0) 72 8080244
 F: +31 (0) 72 8080142

Switzerland 
 Via Lisano 3,
 (CP.683)
 CH-6900
 T: +41 (0) 91 604 55 22
 F: +41 (0) 91 605 17 85

Germany 
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880