



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

LY201 : Lymphoma tissue array, 20 cases/ 20 cores

Microarray Panel	Lymphoma tissue microarray, containing 11 cases of B cell lymphoma, 3 T cell lymphoma, 2 Hodgkin's lymphoma, plus 4 cancer adjacent tissue, single core per case		
Cores	20	LY201 H&E 	
Cases	20		
Row number	4		
Column number	5		
Core Diameter (mm)	2		
Thickness (µm)	5		
Tissue Array Type	FFPE		
Species	Human		
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.		
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "*" means invalid core; "-" means no applicable or negative in IHC markers.		

	1	2	3	4	5
A	Sma	Par	Thy	Ton	Spl
B	Lym	Lym	Lym	Lym	Sto
C	Par	Lym	Lym	Spl	Lym
D	Lym	App	Ton	Spl	Lym

Legend: App - Appendix, Lym - Lymph node, Par - Parotid gland, Sma - Small intestine, Spl - Spleen, Sto - Stomach, Thy - Thyroid gland, Ton - Tonsil
● - Malignant tumor, ● - NAT

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	40	M	Small intestine	B cell lymphoma				malignant	
A2	2	64	M	Parotid gland	B cell lymphoma				malignant	
A3	3	42	F	Thyroid gland	B cell lymphoma				malignant	
A4	4	28	M	Tonsil	B cell lymphoma				malignant	
A5	5	45	M	Spleen	B cell lymphoma				malignant	
B1	6	52	M	Lymph node	B cell lymphoma of right groin				malignant	
B2	7	40	F	Lymph node	Large B cell lymphoma of right oter				malignant	
B3	8	62	F	Lymph node	B cell lymphoma of left armpit				malignant	
B4	9	62	M	Lymph node	B cell lymphoma of right groin				malignant	
B5	10	56	F	Stomach	Mucosa-associated B-cell lymphoma				malignant	
C1	11	51	F	Parotid gland	Mucosa-associated B-cell lymphoma				malignant	
C2	12	63	F	Lymph node	T cell lymphoma of right abdominal part				malignant	
C3	13	55	M	Lymph node	T cell lymphoma of right oter				malignant	
C4	14	43	F	Spleen	T cell lymphoma				malignant	
C5	15	65	F	Lymph node	Sclerotic type Hodgkin's lymphoma of left groin				malignant	
D1	16	4	M	Lymph node	Hodgkin's lymphoma of right neck				malignant	
D2	17	15	F	Appendix	Chronic appendicitis				NAT	

D3	18	35	M	Tonsil	Cancer adjacent normal tonsil tissue				NAT	
D4	19	9	F	Spleen	Cancer adjacent normal spleen tissue				NAT	
D5	20	56	M	Lymph node	Cancer adjacent normal lymph node tissue				NAT	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
 184 Park Drive, Milton Park
 Abingdon OX14 4SE.
 T: +44 (0) 1235 828 200
 F: +44 (0) 1235 820 482

North America 
 1035 Cambridge Street,
 Cambridge, MA 02141.
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

Europe 
 Berenkoog 41,
 1822 BH Alkmaar,
 Netherlands
 T: +31 (0) 72 8080244
 F: +31 (0) 72 8080142

Switzerland 
 Via Lisano 3,
 (CP.683)
 CH-6900
 T: +41 (0) 91 604 55 22
 F: +41 (0) 91 605 17 85

Germany 
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880