



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Balb/C Mouse Total RNA, Panel of 15 Tissues**Catalog #:** MR-015-BLC**Quantity:** 15x0.025mg**Storage Conditions:** Store at -80°C. It is good for one year of the date of purchase if stored properly.**Applications:** RNA is ideal for Northern blotting, ribonuclease protection assay, SI nuclease assay, RT-PCR analysis, rapid amplification of cDNA ends (RACE) and purification of mRNA for library construction.**Quality Control:** The integrity of each RNA sample as indicated by intact 28S and 18S ribosomal RNA is verified by denaturing agarose gel electrophoresis. The purity of RNA is assessed by spectrophotometer (A260/A280 >1.9).**Description:** Total RNA is isolated from 15 mouse tissues (Balb/C strain) selected by the customer using a modified guanidine thiocyanate technique. RNA samples are treated with RNase-free DNase to remove residual DNA, suspended in water, precisely quantified, and stored at -80°C.**List of Tissues Available:**

- | | |
|-------------------------------------|--------------------------------|
| 1) Bladder | 21) Ovary |
| 2) Brain | 22) Pancreas |
| 3) Cecum | 23) Penis |
| 4) Colon | 24) Rectum |
| 5) Diaphragm | 25) Salivary Gland, Mandibular |
| 6) Epididymis, whole | 26) Salivary Gland, Parotid |
| 7) Esophagus | 27) Seminal Vesicles |
| 8) Eye | 28) Skeletal Muscles |
| 9) Heart, left and right atrium | 29) Skin |
| 10) Heart, left and right ventricle | 30) Spinal Cord |
| 11) Heart, whole | 31) Spleen |
| 12) Intestine, Duodenum | 32) Stomach |
| 13) Intestine, Ileum | 33) Testis |
| 14) Intestine, Jejunum | 34) Thymus |
| 15) Intestine, whole | 35) Tongue |
| 16) Kidney | 36) Uterus, non-pregnant |
| 17) Lacrimal Gland | 37) Vagina |
| 18) Liver | |
| 19) Lung | |
| 20) Mammary Gland, non-pregnant | |

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85