



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ST481 : Stomach cancer tissue array with corresponding (self-matching) normal adjacent tissues (12 cases/48 cores), with stage and grade info

Microarray Panel	Stomach carcinoma tissue microarray, containing 12 cases of stomach tumor with their corresponding normal tissues, quadruplicate cores per case	
Cores	48	ST481 H&E 
Cases	12	
Row number	6	
Column number	8	
Core Diameter (mm)	1.5	
Thickness (µm)	5	
Quality Control	Anti-Actin confirmed	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "*" means invalid core; "-" means no applicable or negative in IHC markers.	

	1	2	3	4	5	6	7	8
A	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
B	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
C	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
D	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
E	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
F	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto

Legend: Sto - Stomach
● - Adjacent tissue 1.5 cm away from tumor, ● - Malignant tumor

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	60	M	Stomach	Adenocarcinoma (sparse)	T3N1M0	1		malignant	
A2	2	60	M	Stomach	Adenocarcinoma	T3N1M0	1		malignant	
A3	3	60	M	Stomach	Cancer adjacent normal gastric mucosa tissue	-	-		adjacent	
A4	4	60	M	Stomach	Cancer adjacent normal gastric mucosa tissue	-	-		adjacent	
A5	5	63	M	Stomach	Adenocarcinoma	T3N0M0	1		malignant	
A6	6	63	M	Stomach	Adenocarcinoma	T3N0M0	1		malignant	
A7	7	63	M	Stomach	Cancer adjacent normal gastric tissue	-	-		adjacent	
A8	8	63	M	Stomach	Cancer adjacent normal gastric tissue	-	-		adjacent	
B1	9	63	F	Stomach	Adenocarcinoma	T3N1M0	1 -- 2		malignant	
B2	10	63	F	Stomach	Adenocarcinoma	T3N1M0	1 -- 2		malignant	
B3	11	63	F	Stomach	Cancer adjacent normal gastric tissue (fibrous tissue)	-	-		adjacent	
B4	12	63	F	Stomach	Cancer adjacent normal gastric mucosa tissue	-	-		adjacent	
B5	13	58	M	Stomach	Adenocarcinoma	T3N1M0	1		malignant	
B6	14	58	M	Stomach	Adenocarcinoma	T3N1M0	1		malignant	
B7	15	58	M	Stomach	Cancer adjacent normal gastric tissue	-	-		adjacent	
B8	16	58	M	Stomach	Cancer adjacent normal gastric tissue	-	-		adjacent	
C1	17	73	M	Stomach	Adenocarcinoma	T3N1M0	2		malignant	
C2	18	73	M	Stomach	Adenocarcinoma	T3N1M0	2		malignant	
C3	19	73	M	Stomach	Cancer adjacent normal gastric mucosa tissue (chronic inflammation)	-	-		adjacent	
C4	20	73	M	Stomach	Cancer adjacent normal gastric mucosa tissue (chronic inflammation)	-	-		adjacent	
C5	21	67	M	Stomach	Adenocarcinoma	T3N1M0	2		malignant	

C6	22	67	M	Stomach	Adenocarcinoma (sparse)	T3N1M0	-		malignant	
C7	23	67	M	Stomach	Cancer adjacent normal gastric tissue (fibrous tissue)	-	-		adjacent	
C8	24	67	M	Stomach	Cancer adjacent normal gastric tissue (fibrous tissue)	-	-		adjacent	
D1	25	48	M	Stomach	Papillary adenocarcinoma	T2N1M0	2		malignant	
D2	26	48	M	Stomach	Papillary adenocarcinoma	T2N1M0	2		malignant	
D3	27	48	M	Stomach	Cancer adjacent normal gastric mucosa tissue (chronic inflammation)	-	-		adjacent	
D4	28	48	M	Stomach	Cancer adjacent normal gastric mucosa tissue (chronic inflammation)	-	-		adjacent	
D5	29	45	M	Stomach	Adenocarcinoma	T3N0M0	3		malignant	
D6	30	45	M	Stomach	Adenocarcinoma	T3N0M0	3		malignant	
D7	31	45	M	Stomach	Cancer adjacent normal gastric tissue	-	-		adjacent	
D8	32	45	M	Stomach	Cancer adjacent normal gastric tissue	-	-		adjacent	
E1	33	60	M	Stomach	Adenocarcinoma	T3N3M0	3		malignant	
E2	34	60	M	Stomach	Adenocarcinoma	T3N3M0	3		malignant	
E3	35	60	M	Stomach	Cancer adjacent normal gastric mucosa tissue	-	-		adjacent	
E4	36	60	M	Stomach	Cancer adjacent normal gastric tissue	-	-		adjacent	
E5	37	54	M	Stomach	Adenocarcinoma	T3N2M0	3		malignant	
E6	38	54	M	Stomach	Adenocarcinoma (sparse)	T3N2M0	3		malignant	
E7	39	54	M	Stomach	Cancer adjacent normal gastric tissue (sparse internal gland and smooth muscle)	-	-		adjacent	
E8	40	54	M	Stomach	Cancer adjacent normal gastric tissue	-	-		adjacent	
F1	41	37	F	Stomach	Adenocarcinoma	T3N0M0	3		malignant	
F2	42	37	F	Stomach	Adenocarcinoma	T3N0M0	3		malignant	
F3	43	37	F	Stomach	Cancer adjacent normal gastric tissue (fibrous tissue)	-	-		adjacent	
F4	44	37	F	Stomach	Cancer adjacent normal gastric mucosa tissue (chronic inflammation)	-	-		adjacent	
F5	45	44	M	Stomach	Adenocarcinoma	T3N1M0	3		malignant	
F6	46	44	M	Stomach	Adenocarcinoma	T3N1M0	3		malignant	
F7	47	44	M	Stomach	Cancer adjacent normal gastric mucosa tissue (chronic inflammation)	-	-		adjacent	
F8	48	44	M	Stomach	Cancer adjacent normal gastric mucosa tissue (sparse)	-	-		adjacent	

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880