



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

T009 : Multiple organ cancer test tissue array, with normal tissue, including TNM and pathology grade, 6 cases/24 cores

Microarray Panel	Multiple organ carcinoma test tissue microarray, containing 4 types of organ cancer (stomach, esophagus, small intestine, colon), 2 normal tissue, quadruple cores per case		
Cores	24		
Cases	6		
Row number	3		
Column number	8		
Core Diameter (mm)	1.5		
Thickness (µm)	5		
Quality Control	Anti-Actin confirmed		
Tissue Array Type	FFPE		
Species	Human		
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.		
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet,"*" means invalid core; "-" means no applicable or negative in IHC markers.		

	1	2	3	4	5	6	7	8
A	Sto	Sto	Eso	Eso	Sto	Sto	Eso	Eso
B	Ile	Ile	Col	Col	Ile	Ile	Col	Col
C	Sto	Sto	Col	Col	Sto	Sto	Col	Col

Legend: Col - Colon, Eso - Esophagus, Ile - Ileum, Sto - Stomach

● - Malignant tumor, ● - Malignant tumor (stage I), ● - Malignant tumor (stage II), ● - Malignant tumor (stage III), ● - Normal tissue

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	42	M	Stomach	Neuroendocrine carcinoma	T2N0M0	-	I	malignant	
A2	2	42	M	Stomach	Neuroendocrine carcinoma	T2N0M0	-	I	malignant	
A3	3	47	M	Esophagus	Squamous cell carcinoma	T3N1M0	2	III	malignant	
A4	4	47	M	Esophagus	Squamous cell carcinoma	T3N1M0	2	III	malignant	
A5	5	42	M	Stomach	Neuroendocrine carcinoma	T2N0M0	-	I	malignant	
A6	6	42	M	Stomach	Neuroendocrine carcinoma	T2N0M0	-	I	malignant	
A7	7	47	M	Esophagus	Squamous cell carcinoma	T3N1M0	2	III	malignant	
A8	8	47	M	Esophagus	Squamous cell carcinoma	T3N1M0	2	III	malignant	
B1	9	63	M	Ileum	Adenocarcinoma	T3N0M0	2	II	malignant	
B2	10	63	M	Ileum	Adenocarcinoma	T3N0M0	2	II	malignant	
B3	11	57	M	Colon	Adenocarcinoma	T3N0M0	3	II	malignant	
B4	12	57	M	Colon	Adenocarcinoma	T3N0M0	3	II	malignant	
B5	13	63	M	Ileum	Adenocarcinoma	T3N0M0	2	II	malignant	
B6	14	63	M	Ileum	Adenocarcinoma	T3N0M0	2	II	malignant	
B7	15	57	M	Colon	Adenocarcinoma	T3N0M0	3	II	malignant	
B8	16	57	M	Colon	Adenocarcinoma	T3N0M0	3	II	malignant	
C1	17	35	M	Stomach	Normal stomach tissue	-	-	-	normal	

C2	18	35	M	Stomach	Normal stomach tissue	-	-	-	normal	
C3	19	35	M	Colon	Normal colon tissue	-	-	-	normal	
C4	20	35	M	Colon	Normal colon tissue	-	-	-	normal	
C5	21	35	M	Stomach	Normal stomach tissue	-	-	-	normal	
C6	22	35	M	Stomach	Normal stomach tissue	-	-	-	normal	
C7	23	35	M	Colon	Normal colon tissue	-	-	-	normal	
C8	24	35	M	Colon	Normal colon tissue	-	-	-	normal	
-	0	58	M	Skin	Malignant melanoma (tissue marker)		-		Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
 184 Park Drive, Milton Park
 Abingdon OX14 4SE.
 T: +44 (0) 1235 828 200
 F: +44 (0) 1235 820 482

North America 
 1035 Cambridge Street,
 Cambridge, MA 02141.
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

Europe 
 Berenkoog 41,
 1822 BH Alkmaar,
 Netherlands
 T: +31 (0) 72 8080244
 F: +31 (0) 72 8080142

Switzerland 
 Via Lisano 3,
 (CP.683)
 CH-6900
 T: +41 (0) 91 604 55 22
 F: +41 (0) 91 605 17 85

Germany 
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880