



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product Number	OOIB00013
Product Page	www.avivasysbio.com/normal-llama-serum-ooib00013.html
Name	Normal Llama Serum (OOIB00013)
Host	LLama
Concentration	60 mg/ml
Product Format	Lyophilized
Reconstitution and Storage	Reconstitution: Rehydrate with 2.0 ml of deionized water. Swirl gentle and let stand for up to 2 hours at 18-25C. Centrifuge reconstituted serum to remove any precipitates. Storage: Store lyophilized material at 2-8C. For long term storage after reconstitution, prepare small aliquots and store at -20C. For storage at 2-8C, add a preservative to prevent growth of bacteria. Shelf Life: Expires 1 year from date of receipt. Prepare working dilution prior to use and then discard.
Datasheets/Manuals	Printable datasheet for Normal Llama Serum (OOIB00013)
Additional Information	Normal Llama Serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.
Application Info	This antibody is suitable for all immunoassay applications. The optimal working dilution should be determined by the investigator.
Purification	Lipid extracted and dialyzed against 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Application	ELISA, CHIP, FC, IHC, IP, MS, ARRAY, DB

AVIVA SYSTEMS BIOLOGY manufactures and sells quality antibody products covering genome wide proteins.

This product is for Research Use Only. Not for diagnostic, human, or veterinary use.

Optimal conditions of its use should be determined by end users.

AVIVA SYSTEMS BIOLOGY

6370 Nancy Ridge Dr., Suite 104, San Diego, CA 92121 USA | Tel: (858)552-6979 | info@avivasysbio.com