



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

<b>Product Number</b>	OOMA00058
<b>Product Page</b>	<a href="http://www.avivasysbio.com/troponin-i-cardiac-98-pure-ooma00058.html">www.avivasysbio.com/troponin-i-cardiac-98-pure-ooma00058.html</a>
<b>Name</b>	Troponin I-Cardiac >98% Pure (OOMA00058)
<b>Concentration</b>	1.55mg/ml (OD280nm)
<b>Source</b>	Human cardiac tissue
<b>Product Format</b>	Purified, Liquid
<b>Reconstitution and Storage</b>	Upon receipt, store at -20C. Avoid multiple freeze/thaw cycles.
<b>Datasheets/Manuals</b>	Printable datasheet for <a href="#">OOMA00058</a>
<b>Additional Information</b>	<b>Buffer:</b> 8M Urea, 50mM Tris-HCl, pH 8.0 containing 15mM beta-mercaptoethanol, 1mM EDTA
<b>Application Info</b>	Specific methodologies have not been tested using this product.
<b>Warning</b>	All human source materials have tested negative for HIV1, HIV2, HCV antibodies and HBsAg, HIV-1 antigen, and Syphilis. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.
<b>Purification</b>	Purified
<b>Predicted Species Reactivity</b>	Human

AVIVA SYSTEMS BIOLOGY manufactures and sells quality antibody products covering genome wide proteins.

This product is for Research Use Only. Not for diagnostic, human, or veterinary use.  
Optimal conditions of its use should be determined by end users.