



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

<b>Product Number</b>	OORA00223
<b>Product Page</b>	<a href="http://www.avivasysbio.com/protein-molecular-weight-marker-oora00223.html">www.avivasysbio.com/protein-molecular-weight-marker-oora00223.html</a>
<b>Name</b>	Protein Molecular Weight Marker (OORA00223)
<b>Conjugation</b>	Unconjugated
<b>Concentration</b>	1X
<b>Product Format</b>	Liquid (in 1x Loading Buffer)
<b>Description of Target</b>	This protein molecular weight marker is a cocktail containing seven purified proteins that range in size from 14 to 120 kDa and thus span a molecular weight range common to most proteins and protein subunits. The protein mixture has been formulated to yield well-defined bands after SDS-polyacrylamide gel electrophoresis and staining with Coomassie blue or other common staining techniques. The intensity of the 67 kDa bovine albumin band has been increased to serve as a reference. The product is supplied in 1X loading buffer and is ready-to-use.
<b>Reconstitution and Storage</b>	Store vial at -20C. Expiration date is six (6) months from date of opening container.
<b>Datasheets/Manuals</b>	Printable datasheet for <a href="#">OORA00223</a>
<b>Additional Information</b>	This product is for research use only and is not intended for therapeutic or diagnostic applications.
<b>Application Info</b>	This product is designed for use in accurate estimation of protein molecular weight in SDS-polyacrylamide gel electrophoresis and western blotting. Thaw the cocktail at room temperature. Gently invert vial several times to ensure that the solution is homogeneous and that any precipitated material is re-dissolved. A loading volume of 5-10 ul is appropriate for most gel formats. These are unstained markers optimized for use with Coomassie blue staining, but the loading volume can be adjusted for compatibility with other staining methods (i.e. silver staining). All assays should be optimized by the user.
<b>Application</b>	WB

AVIVA SYSTEMS BIOLOGY manufactures and sells quality antibody products covering genome wide proteins.

This product is for Research Use Only. Not for diagnostic, human, or veterinary use.  
Optimal conditions of its use should be determined by end users.

AVIVA SYSTEMS BIOLOGY

6370 Nancy Ridge Dr., Suite 104, San Diego, CA 92121 USA | Tel: (858)552-6979 | [info@avivasysbio.com](mailto:info@avivasysbio.com)