



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

<b>Product Number</b>	OORA00599
<b>Product Page</b>	<a href="http://www.avivasysbio.com/rabbit-b-non-pooled-complement-oora00599.html">www.avivasysbio.com/rabbit-b-non-pooled-complement-oora00599.html</a>
<b>Name</b>	Rabbit "B" (non-pooled) COMPLEMENT (OORA00599)
<b>Conjugation</b>	Unconjugated
<b>Concentration</b>	85 mg/ml
<b>Product Format</b>	Fresh Frozen
<b>Reconstitution and Storage</b>	Store vial at -70 C prior to opening. NOTE: COMPLEMENT IS A TEMPERATURE SENSITIVE PRODUCT. IMPROPER STORAGE WILL INACTIVATE COMPLEMENT ACTIVITY. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of opening.
<b>Datasheets/Manuals</b>	Printable datasheet for <a href="#">OORA00599</a>
<b>Additional Information</b>	<b>Buffer:</b> None
<b>Additional Information</b>	<b>Label:</b> Unconjugated <b>Preservative:</b> None <b>Stabilizer:</b> None <b>Strain:</b> New Zealand White
<b>Application Info</b>	pH: normalImmunoelectrophoresis: normalHemoglobin: normalIgG Concentration: normalImmunoelectrophoresis: normal
<b>Predicted Species Reactivity</b>	Rabbit

AVIVA SYSTEMS BIOLOGY manufactures and sells quality antibody products covering genome wide proteins.

This product is for Research Use Only. Not for diagnostic, human, or veterinary use.  
Optimal conditions of its use should be determined by end users.

AVIVA SYSTEMS BIOLOGY

6370 Nancy Ridge Dr., Suite 104, San Diego, CA 92121 USA | Tel: (858)552-6979 | [info@avivasysbio.com](mailto:info@avivasysbio.com)