

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com



Sterile Human Blood in Sodium Heparin (OORA01098)

Data Sheet

Product Number	OORA01098
Product Page	www.avivasysbio.com/sterile-human-blood-in-sodium-heparin-oora01098.html
Name	Sterile Human Blood in Sodium Heparin (OORA01098)
Conjugation	Unconjugated
Source	Human. Mixed Sex Adult.
Product Format	Liquid
Reconstitution and Storage	This item is perishable and is collected and processed upon receipt of your order. Store product at 4°C prior to opening. Use aseptic technique to maintain sterility when opening product. Be advised that blood is a perishable product and exact shelf may depend on application. This product MAY be stable for up to four (4) weeks if properly stored and handled.
Datasheets/Manuals	Printable datasheet for OORA01098
Additional Information	Anticoagulant: Sodium Heparin
Warning	No test method can provide total assurance that the hepatitis B virus, hepatitis C virus, human immunodeficiency virus, or any other infectious agents are absent. Thus, all blood products, including purified proteins derived from human blood sources, should be handled at Biosafety Level 2 as recommended by the CDC/NIH manual entitled Biosafety in Microbiological and Biomedical Laboratories for potentially infectious human serum, blood specimens or proteins derived from same. Source material for the human blood product supplied to your facility has been tested for the detection of HIV antibody, Hepatitis B surface antigen, antibody to Hepatitis C, HIV 1 antigen(s), antibody to HTLV - I/II, and syphilis by FDA guidelines. All units were found to be non-reactive/negative for these tests. All human blood source material is collected in FDA licensed centers and is tested with FDA approved test kits.
Preparation	This product was asceptically drawn and processed using sterile equipment and sterile techniques. The product was tested on trypticase soy agar for 24 hours, 48 hours, and 72 hours and was found to be negative for bacteria.
Predicted Species Reactivity	Human

AVIVA SYSTEMS BIOLOGY manufactures and sells quality antibody products covering genome wide proteins.

This product is for Research Use Only. Not for diagnostic, human, or veterinary use. Optimal conditions of its use should be determined by end users.

AVIVA SYSTEMS BIOLOGY 6370 Nancy Ridge Dr., Suite 104, San Diego, CA 92121 USA | Tel: (858)552-6979 | info@avivasysbio.com