



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product Number	OOSA09544
Product Page	www.avivasysbio.com/elisa-coating-buffer-accessory-reagent-oosa09544.html
Name	ELISA COATING BUFFER Accessory Reagent (OOSA09544)
Concentration	5X
Product Format	Liquid
Description of Target	5x ELISA COATING BUFFER
Reconstitution and Storage	2°C to 8°C
Datasheets/Manuals	Printable datasheet for ELISA COATING BUFFER Accessory Reagent (OOSA09544)
Additional Information	Instructions for Use: 1) Dilute the required volume of OOSA09544 1:5 with deionized water and mix for 5-15 minutes to prepare the working solutions. 2) Add antibody or antigen to the working solution and mix for 5-15 minutes. 3) Coat plate (suggested volume 100-200 μ L per well) and incubate at room temperature for 6-24 hours, protected from light and covered to avoid evaporation.
Application	ELISA

AVIVA SYSTEMS BIOLOGY manufactures and sells quality antibody products covering genome wide proteins.

This product is for Research Use Only. Not for diagnostic, human, or veterinary use.
Optimal conditions of its use should be determined by end users.