



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product Number	OOSA09575
Product Page	<a href="http://www.avivasysbio.com/ovalbumin-hrp-accessory-reagent-oosa09575.html">www.avivasysbio.com/ovalbumin-hrp-accessory-reagent-oosa09575.html</a>
Name	OVALBUMIN:HRP Accessory Reagent (OOSA09575)
Protein Size (# AA)	386 amino acids
Conjugation	HRP
NCBI Gene Id	396098
Concentration	1 mg/ml
Alias Symbols	B6.2, Bu-1, Bu-1b
Product Format	Liquid. Phosphate buffered saline
Description of Target	OVALBUMIN:HRP
Reconstitution and Storage	2°C to 8°C -20°C
Datasheets/Manuals	Printable datasheet for <a href="#">OVALBUMIN:HRP Accessory Reagent (OOSA09575)</a>
Additional Information	<b>Approx Protein Conc:</b> 1.0mg/ml
Application Info	<b>ELISA:</b> 1/50 - 1/150
Additional Information	<b>Intended Use:</b> This product is a preparation of HRP conjugated ovalbumin, which may be used for the measurement of ovalbumin specific IgE in immunoassay systems, without species restriction. <b>Preservative:</b> HRP stabiliser (BUF052A) 0.01% Thiomersal
Protein Accession #	<a href="#">NP_990483.1</a>
Gene Symbol	<a href="#">SERPINB14</a>
Application	ELISA

AVIVA SYSTEMS BIOLOGY manufactures and sells quality antibody products covering genome wide proteins.

This product is for Research Use Only. Not for diagnostic, human, or veterinary use.  
Optimal conditions of its use should be determined by end users.