



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product Number	OOSB00005
Product Page	www.avivasysbio.com/gelatin-subbed-slides-oosb00005.html
Name	Gelatin Subbed Slides (OOSB00005)
Product Format	Premium slides subbed with Gelatin
Reconstitution and Storage	Store at -80C
Storage	Gelatin subbed slides are supplied as 72 slides per box. The slides have been carefully cleaned and subbed with a solution of gelatin and chromium potassium sulfate. The slides need to be stored at -70°C temperature until ready to use, at that time they can be thawed.
Datasheets/Manuals	Printable datasheet for OOSB00005
Additional Information	Description: SouthernBiotech's Gelatin Subbed Slides have been carefully cleaned and subbed with a solution of 0.1% gelatin and chromium potassium sulfate. The coated slides measure 25 mm x 75 mm with a frosted end for marking. For optimal adhesion, we recommend sections less than 100 microns in thickness.
Application Info	Mounting sections
Warning	This material contains Chromium (III) Potassium Sulfate Dodecahydrate. Avoid contact with skin, eyes, or clothing. May cause allergic skin reaction. If skin or eye contact occurs, wash with copious amounts of water. If ingested or inhaled, contact a physician immediately.

AVIVA SYSTEMS BIOLOGY manufactures and sells quality antibody products covering genome wide proteins.

This product is for Research Use Only. Not for diagnostic, human, or veterinary use.
Optimal conditions of its use should be determined by end users.