



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

<b>Product Number</b>	OOSB00006
<b>Product Page</b>	<a href="http://www.avivasysbio.com/fluoromount-g-oosb00006.html">www.avivasysbio.com/fluoromount-g-oosb00006.html</a>
<b>Name</b>	Fluoromount-G (OOSB00006)
<b>Product Format</b>	<b>Liquid</b> 25 mL containing sodium azide as a preservative.
<b>Reconstitution and Storage</b>	Store at 22-25°C
<b>Datasheets/Manuals</b>	Printable datasheet for <a href="#">OOSB00006</a>
<b>Additional Information</b>	<p><b>Instruction For Use:</b></p> <ol style="list-style-type: none"> <li>1. Prepare samples and perform staining according to established procedures.</li> <li>2. Following aqueous final step, remove excess liquid around the sample and add 1 drop (~50 uL) of Fluoromount-G for an area of 22x22 mm</li> <li>3. Mount coverslip and remove excess mounting medium.</li> <li>4. Allow to air dry for 5 minutes prior to analysis.</li> </ol>
<b>Application Info</b>	Slide mounting
<b>Warning</b>	Reagents contain sodium azide which is very toxic if ingested or inhaled. Avoid contact with skin, eyes, or clothing. Sodium azide yields toxic hydrazoic acid under acidic conditions. Dilute azide-containing compounds in running water before discarding to avoid accumulation of potentially explosive deposits in lead or copper plumbing.
<b>Application</b>	IHC-F, IHC-P, IHC, ICC, PLA

AVIVA SYSTEMS BIOLOGY manufactures and sells quality antibody products covering genome wide proteins.

This product is for Research Use Only. Not for diagnostic, human, or veterinary use.  
Optimal conditions of its use should be determined by end users.