



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product Number	OPPA02023
Product Page	www.avivasysbio.com/dengue-polyvalent-protein-oppa02023.html
Name	Dengue Polyvalent Protein (OPPA02023)
Protein Size (# AA)	Recombinant amino acids
Purity	Protein is >95% pure as determined by 12% PAGE (coomassie staining).
Source	E. coli
Alias Symbols	Dengue Polyvalent, Polyvalent Dengue Antigen Recombinant
Product Format	Liquid. Supplied in Phosphate buffered saline, pH-7.4 and 0.05% NaN ₃ .
Description of Target	Caused by one of four closely related virus serotypes of the genus Flavivirus, family Flaviviridae, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholino antisense oligos have shown specific activity against Dengue virus.
Reconstitution and Storage	Polyvalent Dengue although stable at 4C for 1 week, should be stored below -18C. Please prevent freeze thaw cycles.
Datasheets/Manuals	Printable datasheet for OPPA02023
Purification	Purified by proprietary chromatographic technique.
Gene Symbol	Dengue Polyvalent
Image 1	

AVIVA SYSTEMS BIOLOGY manufactures and sells quality antibody products covering genome wide proteins.

This product is for Research Use Only. Not for diagnostic, human, or veterinary use.
Optimal conditions of its use should be determined by end users.