



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Rabbit IgG-heavy and light chain Antibody

F(ab')₂ Donkey Polyclonal Conjugate HRP
Antigen Affinity Purified
Catalog No. A120-116P
Lot No. A120-116P-3



APPLICATIONS	WB, IHC, ICC, ELISA								
REACTIVITY TESTED	Rabbit. May cross react with other species.								
ISOTYPE	IgG								
AMOUNT	1 ml at 0.5 mg/ml								
STORAGE/SHELF LIFE	2 - 8° C / 1 year from date of receipt								
PHYSICAL STATE	Liquid								
BUFFER	Phosphate Buffered Saline (PBS) containing 0.2% BSA and 0.1% Pro-Clean 400								
ORIGIN	USA								
PRODUCTION PROCEDURES	<p>The antibody to rabbit IgG was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab')₂ fragments were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to horseradish peroxidase (HRP).</p> <p>Antibody concentration was determined by extinction coefficient: absorbance at 280 nm of 1.4 equals 1.0 mg of IgG.</p> <p>By immunoelectrophoresis and ELISA this antibody reacts specifically with rabbit IgG and with light chains common to other rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with IgG from other species.</p>								
APPLICATIONS	<p>Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.</p> <table><tr><td>Western Blot</td><td>1:5,000 - 1:50,000</td></tr><tr><td>Immunohistochemistry</td><td>1:200 - 1:5000</td></tr><tr><td>Immunocytochemistry</td><td>1:200 - 1:5000</td></tr><tr><td>ELISA</td><td>1:10,000 – 1:100,000</td></tr></table>	Western Blot	1:5,000 - 1:50,000	Immunohistochemistry	1:200 - 1:5000	Immunocytochemistry	1:200 - 1:5000	ELISA	1:10,000 – 1:100,000
Western Blot	1:5,000 - 1:50,000								
Immunohistochemistry	1:200 - 1:5000								
Immunocytochemistry	1:200 - 1:5000								
ELISA	1:10,000 – 1:100,000								
APPLICATION NOTES	Not all listed applications have been specifically tested by our laboratory.								
ADDITIONAL INFO	<p>http://www.bethyl.com/product/A120-116P</p> <p>Use the link above to view SDS, a current list of citations, and other product specific information.</p>								