



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien


T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Detection of human CSN3 by western blot. Samples: Whole cell lysate (50 μ g) from HeLa and HEK293T cells prepared using NETN lysis buffer. **Antibodies:** Affinity purified rabbit anti-CSN3 antibody A300-029A (lot A300-029A-1) used for WB at 0.05 μ g/ml. **Detection:** Chemiluminescence with an exposure time of 3 minutes.



Detection of human CSN3 by western blot. Samples: Whole cell lysate (50 μ g – E; 10 μ g – T) from mock transfected (E) or CSN3 transfected (T) HEKHEK293T cells. **Antibody:** Affinity purified rabbit anti-CSN3 A300-029A used at the indicated concentration. **Detection:** Chemiluminescence with a 10 second exposure.