



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien


T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Detection of human FKBP-135 by western blot and immunoprecipitation. *Samples:* Nuclear extract (NE; 1 mg for IP; 25 μ g for western blot) from HeLa cells. *Antibody:* Affinity purified rabbit anti-FKBP-135 antibodies BL1008 (Cat. No. A300-209A), BL1009 (Cat. No. A300-210A) and BL1010 (Cat. No. A300-211A-T) used at 2 μ g/mg extract for IP. FKBP-135 in NE and immunoprecipitated FKBP-135 were detected using BL1008 (Cat. No. A300-209A) at 0.2 μ g/ml. *Detection:* Chemiluminescence with an exposure time of 10 minutes.