



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien


T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Detection of human CtIP by immunoprecipitation.

Samples: Whole cell lysate (1 mg for IP; 1/4 of IP loaded) from MCF-7 cells. *Antibody:* Affinity purified rabbit anti-CtIP antibody BL1913 (Cat. No. A300-487A-T) used at 3 $\mu\text{g}/\text{mg}$ lysate for IP. CtIP was also immunoprecipitated with rabbit anti-CtIP antibody BL1914 (Cat. No. A300-488A) using 3 $\mu\text{g}/\text{mg}$ lysate. immunoprecipitatesd CtIP was blotted using BL1914 at 1 $\mu\text{g}/\text{ml}$. *Detection:* Chemiluminescence with an exposure time of 30 seconds.