



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien


T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Detection of human CoCoA by western blot of immunoprecipitates. *Samples:* Whole cell lysate (50 μ g for WB; 1 mg for IP, 20% of IP loaded) from MCF-7 cells.

Antibodies: Affinity purified rabbit anti-CoCoA antibody BL2455 (Cat. No. A300-817A) used for IP at 3 μ g/mg lysate (B). CoCoA was also immunoprecipitated by rabbit anti-CoCoA antibody BL2456 (Cat. No. A300-818A), which recognizes a downstream epitope. For blotting immunoprecipitated CoCoA, BL2456 was used at 1 μ g/ml.

Detection: Chemiluminescence with an exposure time of 30 seconds.