



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Normal Goat Serum

Goat

Catalog No. A50-001S-100ML

AMOUNT	100 ml
STORAGE/SHELF LIFE	2 – 8°C / 1 year from date of receipt
PHYSICAL STATE	Liquid
BUFFER	Serum containing 0.09% Sodium Azide
ORIGIN	USA
PRODUCTION PROCEDURES	Goat serum is collected from normal, healthy animals.

APPLICATIONS Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.

APPLICATION NOTES Normal goat serum can be used for non-specific blocking of surfaces, tissue, and protein in applications such as IHC, western blot, and ICC. This generally occurs by binding to Fc receptors that could interact with primary or secondary antibodies, thus causing off-target binding, or simply interacting with any other surfaces that may non-specifically bind IgGs.

ADDITIONAL INFO <https://www.fortislife.com/p/A50-001S-100ML>
Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.
Michael Spencer, PhD Date: March 1, 2023