



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# PTEN Recombinant Monoclonal Antibody [BLR270L]

Rabbit Recombinant Monoclonal

Purified RefSeq ID NP\_000305.3  
Catalog No. A700-270CF Uniprot ID P60484  
Lot No. 1

---

**SPECIES REACTIVITY** Human, Mouse, Rat  
**AMOUNT** 100 µl  
**CONCENTRATION** 1000 µg/ml  
**STORAGE/SHELF LIFE** 2 – 8°C / 1 year from date of receipt  
**PHYSICAL STATE** Liquid  
**BUFFER** Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, BSA-Free  
**ISOTYPE** IgG  
**CLONE #** BLR270L  
**ORIGIN** USA  
**PRODUCTION PROCEDURES** Recombinant antibody was purified from cell culture supernatant.

Immunogen was a peptide representing a region between residue 353 and the C-terminus (residue 403) of human Phosphatidylinositol 3,4,5-trisphosphate 3-phosphatase and dual-specificity protein phosphatase PTEN using the numbering given in entry NP\_000305.3 (Gene ID 5728).

**APPLICATIONS** Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.

A700-270CF is the carrier-free version of 2772-1C7, which is qualified for use in Western Blot, Immunoprecipitation, Immunohistochemistry, and/or Immunocytochemistry. The format of A700-270CF is designed for compatibility with the labeling of the antibody such as with fluorochromes, metal isotopes, oligonucleotides, and enzymes. Upon completion of labeling, the user must empirically define the assay dependent concentration for use.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.  
Michael Spencer, PhD Date: June 30, 2023