



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Calretinin, monoclonal

nordicmubio.com/product/calretinin-monoclonal

Catalogue number: **CA130**

Clone	6B3
Isotype	ig
Product Type	Primary Antibodies
Units	1 ml
Host	Mouse
Application	Immunoblotting Immunohistochemistry (frozen) Immunohistochemistry (paraffin)

Background

The antibody is suitable for the detection of Calretinin in neurons and may be used as an aid in differential diagnosis of pleura-mesothelioma. This marker is frequently expressed in epitheloid mesothelioma, while occasionally it turns positive in sarcomatous mesothelioma and carcinoma. It was found to be negative in sarcoma (2,3) and in 90% of adenocarcinoma. Calretinin is a Calcium-binding protein related to Calmodulin and Calbindin. It is found in groups of neurons of the brain cortex, in the retina and in the spinal cord (1). More recent studies show that it is a suitable marker for mesothelioma of the lung (2, 3, 4, 5). Calretinin, no cross-reaction with Calbindin and Calmodulin.

Source

Immunogen: Recombinant human Calretinin 22kDa

Product

Antibody solution in stabilizing phosphate buffer pH 7.3. Contains 0.09 % sodium azide**. The volume is sufficient for 200-500 immunohistochemical tests (100 µl working solution / test). Use appropriate antibody diluent e.g. BIOLOGO Art. No. PU002, if further dilution is required.

Purification Method: Antibody solution in stabilizing phosphate buffer pH 7.3. Contains 0.09 % sodium azide**. The volume is sufficient for 200-500 immunohistochemical tests (100 µl working solution / test). Use appropriate antibody diluent e.g. BIOLOGO Art. No. PU002, if further dilution is required.

Secondary Reagents: We recommend the use of BIOLOGO's Universal Staining System DAB (Art. No. DA005) or AEC (Art. No. AE005).

Specificity

Species Reactivity: Human, Rat

Applications

IHC(C, P), Immunoblotting

Incubation Time: 60 min at RT

Working Concentration: (liquid conc.) 1:40-1:50

Pre-Treatment: Improvement of staining with Unmasking Fluid G (Art. No. DE007) or Unmasking Fluid C (Art. No. DE000)

Positive Control: Cerebral cortex, pleural mesothelioma

Storage

2-8°C

Caution

*These antibodies are intended for in vitro research use only. They must not be used for clinical diagnostics and not for in vivo experiments in humans or animals. ** The preservative sodium azide is known to be poisonous and potentially hazardous to health. It should be handled only by trained staff. Despite of the product's low azide concentration it must be handled with care. Dispose according to regional rules!

References

1. Schwaller B., Buchwald P., Blümcke I., Celio M.R., and Hunziker W. (1994) Charakterization of a polyclonal antiserum against the purified human recombinant calcium binding protein calretinin. *Cell Calcium* 14; 639-648. 2. Gotzos V., Vogt P., and Celio M.R. (1995) Calretinin is a selective marker for malignant pleural mesotheliomas of the epithelial type. *Pathol. Res. Pract.* 192; 137-147. 3. Doglioni C., Tos A.P., Laurino L., Iuzzolino P., Chiarelli C., Celio M.R., and Viale G. (1996) Calretinin: a novel immunohistochemical marker for mesothelioma. *Am. J. Surg. Pathol.* 20; 1037-1046. 4. Zimmermann L., and Schwaller B. (2002) Monoclonal antibodies recognizing epitopes of calretinins: dependence on Ca²⁺-binding status and differences in antigen accessibility in colon cancer cells. *Cell Calcium* 31; 13-25.

Protein Reference(s)

Database Name: UniProt

Accession number: P22676 (CALB2_HUMAN)

Species Accession: Human