



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



---

## Collagen Type III, chicken

---

[nordicmubio.com/product/collagen-type-iii-chicken-2](http://nordicmubio.com/product/collagen-type-iii-chicken-2)

Catalogue number: **CO20331-0.1**

Product Type	Primary Antibodies
Units	0.1 ml
Host	Rabbit
Application	ELISA Immunofluorescence Immunohistochemistry (paraffin) Radioimmunoassay Western Blotting

### Background

Type III collagen is an  $\alpha_1(\text{III})$ -trimer, which forms 67 nm cross-banded fibrils. Typically it can be observed in skin, cartilage and vitreous body. Collagens consist of a family of highly specialized glycoproteins of which at least 16 genetically distinct types are known to date. The basal unit of a collagen molecule consists of a triple-helical structure formed by 3 alpha-chains. Predominant amino acids are glycine, proline and hydroxyproline. Regularly also lysines and hydroxylysines occur, which are responsible for cross-linkage and glycosylation of the protein chains. Different composition of alpha-chains and different glycosylation contribute to the high variability of collagens in different tissues and organs. Chicken collagen type III 100%, chicken collagen type I <5%.

### Source

*Immunogen:* Purified collagen type III from chicken skin

### Product

affinity purified antibody lyophilized from phosphate buffered solution; no BSA and preservative added!

*Purification Method:* affinity purified antibody lyophilized from phosphate buffered solution; no BSA and preservative added!

*Concentration:* app. 1 mg/ml

*Secondary Reagents:* Anti-rabbit IgG-conjugates, e.g. anti-rabbit IgG:FITC (Art. No. FI-1000) or anti-rabbit IgG:DyLight488 (Art. No. DI-1488).

## **Specificity**

*Species Reactivity:* Chicken, cross-reaction with human, bovine, mouse and rat collagen type III <0,1%

## **Applications**

IHC(P), IFA, ELISA, RIA, IB/WB

*Incubation Time:* IHC(P) 60 min at RT or 2-8°C over night

*Working Concentration:* (purified, lyophilized) IFA 1:40, IHC(P) 1:500, ELISA > 1:200

*Pre-Treatment:* After de-waxing the tissue slices they are treated with 0.2% hyaluronidase (app. 300 U/mg e.g. Art. No. HYA02-50) in TBS 15 min at 37°C. Thereafter non-specific binding is blocked by blocking serum or 3% BSA in TBS. For peroxidase systems blocking with 1% peroxide solution in TBS for 30 min at RT is recommended.

*Positive Control:* Chicken skin or liver

## **Storage**

-20°C

## **Caution**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. It may contain hazardous ingredients. Please refer to the Safety Data Sheets (SDS) for additional information and proper handling procedures. Dispose product remainders according to local regulations. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.

## **References**

1. Robert J., Hartmann D.J., Sengel P. (1989) Production of fibronectin and collagen types I and III by chick embryo dermal cells cultured on extracellular matrix substrates. *Int. J. Dev. Biol.* 33, 267-275.
2. Robert J., Mauger A., Sengel P. (1989) Influence of various extracellular matrix components on the behavior of cultured chick embryo dermal cells. *Int. J. Dev. Biol.* 33, 227-237.
3. Duband J.L., Thiery J.P. (1987) Distribution of laminin and collagens during avian neural crest development. *Development* 101, 461-478.
4. Mauger A., Demarchez M., Herbage D., Grimaud J.A., Druguet M., Hartmann D.J. Foidart J.M., Sengel P. (1983) Immunofluorescent localization of collagen types I, III, IV, fibronectin and laminin during morphogenesis of scales and scaleless skin in the chick

**Protein Reference(s)**

*Database Name:* UniProt

*Accession number:* P12105 (CO3A1\_CHICK)

*Species Accession:* Chicken