



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Leukemia Inhibitory Factor, His-GST tag, human recombinant (rHuLIF-His-GST)

Catalog No: 97532
Lot No: XXXXX
Source: *E. coli*
Synonyms: D factor, MLPLI, HILDA, Emfilermin, Leukemia Inhibitory factor, Differentiation-stimulating factor, Melanoma-derived LPL inhibitor

Background

Leukemia inhibitory factor, is a pleiotropic cytokine which is expressed by numerous cells including activated T lymphocytes, monocytes, mast cells and neuronal cells. LIF takes part in the induction of hematopoietic differentiation in normal and myeloid leukemia cells, induction of neuronal cell differentiation, regulator of mesenchymal to epithelial conversion during kidney development, and is a key player in immune tolerance at the maternal-fetal interface.

Description

Leukemia Inhibitory Factor (LIF) produced in *E. coli* is a single, non-glycosylated polypeptide chain containing 415 amino acids (23-202) and having a molecular mass of 47.2 kDa. LIF is fused to a 236 amino acid His-GST tag at N-terminus and purified by proprietary chromatographic techniques.

Physical Appearance

Sterile filtered clear solution.

Formulation

The LIF GST protein solution (1 mg/1 ml) is formulated in 20 mM Tris-HCl buffer (pH 8.0) 1 mM DTT, 50 mM NaCl and 20% glycerol.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.

Purity

Greater than 90% as determined by SDS-PAGE.

Amino Acid Sequence

MHHHHHMSPLILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDV KLTQSMAIIR
YIADKHNMLG GCPKERA EIS MLEGAVLDIR YGVSRIAYSK DFETLKVDFL SKLPEMLKMF EDRLCHKTYL NGDHVTHPDF
MLYDALDVVL YMDPMCLDAF PKLVCFKKRI EAIPQIDKYL KSSKYIAWPL QGWQATFGGG DHPPKSDLVP RGSMSPLPI
TPVNATCAIR HPCHNNLMNQ IRSQLAQLNG SANALFILYY TAQGEPPFNN LDKLCGPNTV DFPPFHANGT EKAKLVELYR
IVVYLGTSLG NITRDQKILN PSALSLHSKL NATADILRGL LSNVLCRLCS KYHVGHVDTV YGPDTSKGDV FQKKKLGCCQL
LGKYKQIIAV LAQAF

Usage

This product is offered by Biomol for research purposes only. Not for diagnostic purposes or human use. It may not be resold or used to manufacture commercial products without written approval of Biomol GmbH.