



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Data Sheet

### ROR2, Fc-Fusion (IgG1), Avi-Tag, Biotin-Labeled

Human, Recombinant, C-terminal Fc-Avi-tag,  
Biotin-Labeled

**Catalog #:** 100046

**Lot #:** 180426      **Conc.:** 1.23 mg/ml

**Description:** Human ROR2, also known as tyrosine-protein kinase transmembrane receptor, neurotrophic tyrosine kinase, and receptor related 2, GenBank Accession No. NM\_004560, a.a. 28-403, with C-terminal Avi-Tag fused to the Fc region of Human IgG1, and expressed in a Hek293 expression system. MW=71 kDa. This protein is enzymatically biotinylated using Avi-Tag™ technology.

**Formulated in:** 8 mM Phosphate, pH 7.4, 110 mM NaCl, 2.2 mM KCl, and 20% glycerol

**Stability:** At least 6 months at -80°C. Avoid freeze/thaw cycles. Storing diluted enzyme is not recommended, if necessary, use carrier protein (BSA 0.1 – 0.5%).

**References:**

1. DeChiara, T.M. *et al.*, *Nature Genet.* 2000; **24**: 271-274.
2. Nomi, M. *et al.*, *Molec. Cell Biol.* 2001; **21**: 8239-8335.


## **Quality Assurance**

**4-20% SDS-PAGE  
Coomassie staining**

**Lane 1:**  
4µg ROR2

**Lane 2:**  
Protein Marker

**MW:**  
71 kDa +glycans  
**Purity:** ≥90%



OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**


Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)



6042 Cornerstone Court West, Suite B  
San Diego, CA 92121  
Tel: 1.858.829.3082  
Fax: 1.858.481.8694  
Email: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

### Gel Filtration Curve



<10% aggregation

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**  
Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)  
Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)