



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Product Information

<b>Description:</b>	Recombinant human KEAP1 (kelch like ECH associated protein 1), full length, encompassing amino acids 2-624. This construct contains an N-terminal His-tag (6xHis). This recombinant protein was affinity purified.
<b>Species:</b>	Human
<b>Construct:</b>	KEAP1 (His-2-624(end))
<b>Concentration:</b>	0.25 mg/ml
<b>Expression System:</b>	<i>E. coli</i>
<b>Purity:</b>	80%
<b>Format:</b>	Aqueous buffer solution
<b>Formulated In:</b>	40 mM Tris-HCl pH 8.0, 110 mM NaCl, 2.2 mM KCl, 20% glycerol
<b>MW:</b>	70 kDa
<b>Genbank Accession:</b>	NM_203500
<b>Stability:</b>	At least 6 months at -80°C.
<b>Storage:</b>	-80°C
<b>Instructions for Use:</b>	Thaw on ice and gently mix prior to use. DO NOT VORTEX. Perform a quick spin before opening. Aliquot into small volumes and flash freeze for long term storage. Avoid multiple freeze/thaw cycles.

## Quality Control Data

### 4-20% SDS-Page Coomassie Staining

