



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com


[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product Information


| | |
|---------------------------|---|
| Construct: | Trop2 (88-274-Fc(IgG1)-Avi)-(Biotin) |
| Label: | This protein is enzymatically biotinylated using Avi-Tag™ technology. Biotinylation confirmed to be ≥90%. |
| Concentration: | 1.02 mg/ml |
| Species: | Human |
| Formulated In: | 8 mM phosphate, pH 7.4, 110 mM NaCl, 2.2 mM KCl, and 20% glycerol |
| Expression System: | HEK293 |
| Format: | Aqueous buffer solution |
| Stability: | At least 6 months at -80°C. Avoid freeze/thaw cycles. |
| Storage: | -80°C |
| Genbank Accession: | NM_002353 |
| MW: | 50 kDa + glycans |
| Glycosylation: | This protein runs at a higher MW by SDS-PAGE due to glycosylation. |
| Purity: | ≥90% |

Quality Control Data

4-20% SDS-Page Coomassie Staining



Biotin-Avidin Pulldown



1. Beads
2. Flow thru
3. Control
4. Standards

* Avidin from beads.