



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Product Information

<b>Description:</b>	Recombinant human DLL1 (delta like canonical Notch ligand 1), encompassing amino acids 18-540. This construct contains a C-terminal Avi-Tag™ followed by a His-tag (6xHis). The recombinant protein was affinity purified.
<b>Background:</b>	DLL1 is involved in cell-to-cell communication and in embryonic development, neurogenesis, and tissue homeostasis. It is a member of the Delta-Serrate-Lag2 (DSL) family of ligands that activate the Notch signaling pathway. DLL1 binds to the Notch receptors on neighboring cells, triggering signals that regulate cell fate determination, proliferation, and differentiation. Dysregulation of DLL1 has been implicated in several diseases.
<b>Species:</b>	Human
<b>Construct:</b>	DLL1 (18-540-Avi-His)
<b>Concentration:</b>	0.12 mg/ml
<b>Expression System:</b>	HEK293
<b>Purity:</b>	≥90%
<b>Format:</b>	Aqueous buffer solution.
<b>Formulated In:</b>	8 mM phosphate, pH 7.4, 110 mM NaCl, 2.2 mM KCl, and 20% glycerol
<b>MW:</b>	60 kDa + glycans
<b>Glycosylation:</b>	This protein runs at a higher MW by SDS-PAGE due to glycosylation.
<b>Genbank Accession:</b>	NM_005618.4
<b>Stability:</b>	At least 6 months at -80°C.
<b>Storage:</b>	-80°C
<b>Instructions for Use:</b>	Thaw on ice and gently mix prior to use. DO NOT VORTEX. Perform a quick spin before opening. Aliquot into small volumes and flash freeze for long term storage. Avoid multiple freeze/thaw cycles.

## Quality Control Data

### 4-20% SDS-PAGE Coomassie Staining

