



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6044 Cornerstone Court West, Suite E
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Sheet

SIRT Assay Buffer

Catalog #: 50090

Lot#: 110727

SIZE: 20 mL

DESCRIPTION: Aqueous buffer optimized for use with BPS Bioscience's Sirtuin 1,2,& 3 (SIRT) assay kits:

- | | | |
|-------------------------------|--------|------------|
| • Fluorogenic SIRT1 Assay Kit | #50081 | 100 assays |
| • Fluorogenic SIRT2 Assay Kit | #50082 | 100 assays |
| • Fluorogenic SIRT3 Assay Kit | #50083 | 100 assays |

APPLICATIONS: Suitable for use with cell extracts, nuclear extracts, immunoprecipitates, and purified SIRT enzymes.

SUPPLIED AS: Tris-buffered solution

PURITY: ≥98% purity by HPLC.

STORAGE: Store at -20°C. Avoid freeze/thaw cycles.

STABILITY: Stable for at least 6 months when stored as recommended.

RELATED PRODUCTS:

HDAC Assay Developer	#50030	6 mL
SIRT Assay Developer	#50089	6 mL
SIRT1 (Sir2) Assay Kit	#50081	100 rxns.
SIRT2 Assay Kit	#50082	100 rxns.
SIRT3 Assay Kit	#50083	100 rxns.
SIRT1 Active Enzyme	#50012	100 µg
SIRT2 Active Enzyme	#50013	100 µg
SIRT3 Active Enzyme	#50014	100 µg
SIRT4 Active Enzyme	#50015	100 µg
SIRT5 Active Enzyme	#50016	100 µg
SIRT6 Active Enzyme	#50017	100 µg
SIRT7 Active Enzyme	#50018	100 µg

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com