



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court W, Ste B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Thaw Medium 5

Catalog #: 60182

Description

This liquid Cell Thawing Medium is optimized for thawing and plating select BPS Bioscience recombinant cell lines.

Formulation

DMEM medium (Hyclone #SH30243.01) supplemented with 10% Calf Serum (CS) (Hyclone # SH30072.03), not FBS, 1% Penicillin/Streptomycin (Hyclone # SV30010.01).

Quality

Sterile filtered (0.2 μ m) and Mycoplasma-free. Each lot is tested for growth of recombinant cells. All components are low endotoxin. CS is ISIA-certified and sourced from the United States of America.

Storage and Stability

Upon receipt, the medium should be stored in the dark at 4°C. Avoid repeated freeze/thaw cycles. Do not exceed 37°C. When stored in the dark at 4°C, the product is stable for six months after first thaw.

Thawing Protocol

For detailed thawing protocol, see cell line data sheet.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: [**info@bpsbioscience.com**](mailto:info@bpsbioscience.com)

Please visit our website at: [**www.bpsbioscience.com**](http://www.bpsbioscience.com)