



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6405 Mira Mesa Blvd Ste 100  
San Diego, CA 92121  
Tel: 1.858.202.1401  
Fax: 1.858.481.8694  
Email: support@bpsbioscience.com

## Growth Medium 3A Catalog #: 60188

### Description

This liquid Growth Medium is a complete medium optimized to culture select BPS Bioscience cell lines.

### Formulation

BPS Growth Media is formulated through BPS Thaw Media combined with antibiotics for maintenance of selective markers in the stable cell line.

- **Growth Medium 3A (BPS #60188):** Thaw Medium 3 (BPS #60186) plus 1 mg/ml of Geneticin (Thermo Fisher, #11811031) and 500 µg/ml of Hygromycin B (Thermo Fisher, #10687010).
  - **Thaw Medium 3 (BPS #60186):** F-12K Medium (Kaighn's Modification of Ham's F-12 Medium) supplemented with 10% FBS (Thermo Fisher, #26140079), 1% Penicillin/Streptomycin (Hyclone SV30010.01).

### Quality

Sterile filtered (0.2 µm) and Mycoplasma-free. Each lot is tested for growth of recombinant cells. All components are low endotoxin. FBS is ISIA-certified and sourced from the United States of America.

### Storage and Stability

Upon receipt, the medium should be stored in the dark at 4°C. Avoid repeated freeze/thaw cycles. Do not exceed 37°C. When stored in the dark at 4°C, the product is stable for six months after first thaw.

### Culture Protocol

For detailed thawing and cell culture protocol, see cell line data sheet.

### Optimized for Use With

Product	Cat. #	Size
PD-L1 / TCR Activator - CHO Recombinant Cell Line	#60536	2 vials
Woodchuck PD-L1 / TCR activator - CHO Recombinant Cell line	#79457	2 vials

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: [support@bpsbioscience.com](mailto:support@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)