



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Data Sheet

OX40, Fc-Fusion, Avi-Tag (Mouse)

Mouse, Recombinant, C-terminal Fc-Avi-Tag
Catalog #: 72546
Lot #: 171004 **Conc.:** 0.53 mg/ml

Formulated in: 8mM Phosphate, 110mM NaCl, 2.2mM KCl, pH7.4, and 20% glycerol.

Stability: At least 6 months at -80°C . Avoid freeze/thaw cycles. Protein may be diluted to $\geq 100 \mu\text{g/ml}$ in PBS + glycerol and stored at -80°C .

Description: Mouse OX40, also known as Tumor Necrosis Factor Receptor Superfamily Member 4, TNFRSF4, and CD134, GenBank Accession No. NM_011659, a.a. 20-211 with C-terminal Avi-Tag fused to the Fc portion of Human IgG1 and expressed in a HEK293 cell expression system. MW = 52 kDa. This protein runs at a higher MW by SDS-PAGE due to glycosylation.

References:

1. Linton, P.J., *et al.*, *J. Exp. Med.* 2003; **197(7)**: 875-883.
2. Arestides, R.S.S., *et al.*, *Eur. J. Immunol.* 2002; **32(10)**: 2874-2880.


Quality Assurance

4-20% SDS-PAGE Coomassie staining

Lane 1:
4 μg OX40(Mouse)

Lane 2:
Protein Marker

MW: 52 kDa
Purity: $\geq 90\%$




OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court W, Ste B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Gel Filtration Curve



<10% Aggregation

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com