



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Data Sheet

### CD111, Fc-Fusion, Avi-Tag, Biotin-Labeled

Human, Recombinant, C-terminal Fc-Avi-Tag


**Catalog #:** 77002**Lot #:** 170831**Conc:** 0.97 mg/ml**Formulated in:** 8 mM Phosphate, pH 7.4, 400 mM NaCl, 2.2 mM KCl, and 20% glycerol.**Stability:** At least 6 months at  $-80^{\circ}\text{C}$ . Avoid freeze/thaw cycles. Storing diluted enzyme is not recommended, if necessary, use carrier protein (BSA 0.1 – 0.5%).**References:**

1. Friedman, G., *et al.*, *Neuro-Oncology*. 2009; **95(2)**; 199-209.
2. Waters, A., *et al.*, *Transitional Oncology*. 2016; **9(5)**: 419-430.

**Description:**

Human CD111, also known as Poliovirus receptor-related 1, netcin-1, herpesvirus entry mediator C, and HVEC. GenBank Accession No. NM\_002855, a.a. 31-355, fused at the C-terminus to the Fc portion of Human IgG1 and with a C-terminal Avi-tag, expressed in a Hek293 cell expression system. This protein is enzymatically biotinylated using Avi-Tag™ technology. MW=65 kDa.

### Quality Assurance

**4-20% SDS-PAGE  
Coomassie staining****Lane 1:**  
4 $\mu\text{g}$  CD111**Lane 2:**  
Protein Marker**MW:**  
65 kDa + glycans**Purity:**  $\geq 90\%$ 


OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** or Fax **1.858.259.5734**Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)



6042 Cornerstone Court West, Suite B  
San Diego, CA 92121  
Tel: 1.858.202.1401  
Fax: 1.858.481.8694  
Email: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

### Gel Filtration Curve



<10% Aggregation

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** or Fax **1.858.259.5734**  
Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)  
Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)