



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Data Sheet

PDPN, Fc-Fusion, Avi-Tag, Biotin-Labeled

Human, Recombinant, C-terminal Fc-Avi-Tag,
Biotin-Labeled

Catalog #: 79342

Lot#: 171026

Conc.: 2.28 mg/ml

Formulated in: 8 mM Phosphate, pH 7.4,
110 mM NaCl, 2.2 mM KCl, and 20%
glycerol.

Stability: At least 6 months at -80°C . *Avoid
freeze/thaw cycles. Protein may be diluted to
 $\geq 100 \mu\text{g/ml}$ in PBS + glycerol and stored at
 -80°C .*

References:

1. Wicki, A., et al., Br. J. Cancer. 2007;
96(1): 1-5.
2. Robertson, N., et al., Genomics.
1995; 23(1):42-50.

Description:

Human PDPN, known as Podoplanin, GenBank Accession No. NM_006474, a.a. 23-131, with C-terminal Avi-Tag fused to the Fc portion of Human IgG1 and expressed in a HEK293 cell expression system. This protein is enzymatically biotinylated using Avi-Tag™ technology. This protein runs at a higher MW by SDS-PAGE due to glycosylation. MW = 40 kDa.

Quality Assurance


4-20% SDS-PAGE Coomassie staining

Lane 1:
4 μg PDPN

Lane 2:
Protein Marker

MW: 40 kDa + glycans

Purity: $\geq 90\%$



OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**


Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court W, Ste B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Gel Filtration Curve



<10% Aggregation

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com