



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Data Sheet

### DLK, GST-Tag

Human, Recombinant, with N-terminal GST-tag

**Catalog #:** 79345

**Lot #:** 171023      **Conc.:** 0.1 mg/ml

**Description:** Recombinant Human DLK, also known as MAP3K12, MUK, and ZPK. Genbank Accession #NM\_006301, a.a. 1-520, with an N-terminal GST-tag, MW = 61 kDa, expressed in Sf9 cells via a baculovirus expression system.

**Formulated in:** 50 mM Tris-HCl, pH 7.5, 50 mM NaCl, 10mM glutathione, 0.1 mM EDTA, 0.25 mM DTT, 0.1 mM PMSF, 25% glycerol.

**Stability:** At least 6 months at -80°C.  
Avoid freeze/thaw cycles.

#### References:

1. Hirai, S., *et al.*, *Oncogene* 2011;**12**:641-650.
2. Ferron, S., *et al.*, *Nature*. 2011;**475**:381-385.

## Quality Assurance

### 10% SDS-PAGE Coomassie staining

**Lane 1:**


DLK

**Lane 2:**

Protein Marker

**MW:** 94 kDa

**Purity:** ≥75%



OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)