



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product Information

Description:	NK Cell Culture Cytokine Cocktail is an optimized 1000x cytokine mixture designed to support primary NK cell culture with improved transduction efficiency, cell viability and expansion after CAR (chimeric-antigen receptor) engineering, in combination with NK Medium, Serum-Free (#82615).
Quality:	Mycoplasma-free. All components are low endotoxin, suitable for cell culture.
Storage:	Upon receipt, the cocktail should be stored at -20°C. After thawing on ice, the cocktail should be aliquoted and stored at -20°C.
Stability:	Do not exceed 37°C. When stored at -20°C, the product is stable for six months. Avoid repeated freeze/thaw cycles.

Directions for Use:

1. Dilute NK Cell Culture Cytokine Cocktail to 1x in NK Medium, Serum-Free. Mix thoroughly, filter if preferred. Set up cultures as desired.