



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court W, Ste B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Sheet

Empty vector control - HEK293 Recombinant Cell line Catalog #: 90334

Product description

Recombinant HEK293 cell transfected with empty expression vector containing the neomycin resistance gene. This cell line serves as a negative control for

- TMEM16A-HEK293 Cell Line Cat. #90230
- TRPC3-HEK293 Cell Line Cat. #90130
- TRPC7-HEK293 Cell Line Cat. #90030
- IKCA1 (KCNN4) - HEK293 Cell Line Cat. #90330
- ERG (Kv11.1) - HEK293 Cell Line Cat. #60619

Format

Each vial contains 1.5×10^6 cells in 1 ml of 10% DMSO.

Storage

Immediately upon receipt, store in liquid nitrogen.

Mycoplasma Testing

The cell line has been screened using the PCR-based Venor™ GeM Mycoplasma Detection kit (Sigma-Aldrich) to confirm the absence of Mycoplasma species.

Culture conditions

Thaw Medium 1 (BPS Cat. #60187): MEM medium (Hyclone #SH30024.01) supplemented with 10% FBS (Invitrogen #26140-079), 1% non-essential amino acids (Hyclone #SH30238.01), 1 mM Na pyruvate (Hyclone #SH30239.01), 1% Penicillin/Streptomycin (Hyclone SV30010.01)

Growth Medium 1B (BPS Cat. #79531): Thaw Medium 1 (BPS Cat. #60187) plus 400 µg/ml of Geneticin (life Technologies #11811031) to ensure the recombinant expression is maintained.

Cells should be grown using at 37°C with 5% CO₂ using Growth Medium 1B. hERG-HEK293 cells should exhibit a typical cell division time of ~24 hours.

To thaw the cells, it is recommended to quickly thaw the frozen cells from liquid nitrogen in a 37°C water-bath and transfer to a tube containing 10 ml of Thaw Medium 1 (**no Geneticin**). Spin down cells, resuspend cells in pre-warmed Thaw Medium 1 (**no Geneticin**), and transfer resuspended cells to T25 flask. Culture at 37°C in a CO₂ incubator overnight. The next day, replace the medium with fresh Thaw Medium 1 (**no Geneticin**) and continue growing culture in

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court W, Ste B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

a CO₂ incubator at 37°C until the cells are ready to be split. Cells should be split before they reach complete confluence. At first passage, switch to Growth Medium 1B (**contains Geneticin**).

To passage the cells, rinse cells with phosphate buffered saline (PBS) and detach cells from culture vessel with 0.05% Trypsin/EDTA. Add Growth Medium 1B and transfer to a tube, spin down cells, resuspend cells, and seed appropriate aliquots of cell suspension into new culture vessels. Subcultivation ration: 1:10 to 1:20 weekly or twice a week.

To freeze cells, rinse cells with phosphate buffered saline (PBS) and detach cells from culture vessel with 0.05% Trypsin/EDTA. Add Growth Medium 1B and transfer to a tube, spin down cells and resuspend in freezing medium (90% FBS, 10% DMSO). Freeze cells at -80°C overnight and transfer to liquid nitrogen the next day.

License Disclosure

Purchase of this cell line grants you with a 10-year license to use this cell line in your immediate laboratory, for research use only. This license does not permit you to share, distribute, sell, sublicense, or otherwise make the cell line available for use to other laboratories, departments, research institutions, hospitals, universities, or biotech companies. The license does not permit use of this cell line in humans or for therapeutic or drug use. The license does not permit modification of the cell line in any way. Inappropriate use or distribution of this cell line will result in revocation of the license and result in an immediate cease of sales and distribution of BPS products to your laboratory. BPS does not warrant the suitability of the cell line for any particular use, and does not accept any liability in connection with the handling or use of the cell line. Modifications of this cell line, transfer to another facility, or commercial use of the cells may require a separate license and additional fees; contact sales@bpsbioscience.com for details. Publications using this cell line should reference BPS Bioscience, Inc., San Diego.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com