



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# Product Information

## Donkey Anti-Guinea Pig IgG (H+L), Whole Antibody, Highly Cross-Adsorbed (min x Bv, Ch, Gt, Hs, Hu, Ms, Rb, Rt, Sh, SHm)

CF® Dye	Ex/Em	Catalog No.	Size
CF®405M	408/452 nm	20376-50uL	50 uL
		20376-500uL	500 uL
		20376-1mg	1 mg (lyophilized)
CF®405S	404/431 nm	20356-50uL	50 uL
		20356-500uL	500 uL
		20356-1mg	1 mg (lyophilized)
CF®488A	490/515 nm	20169-1	50 uL
		20169	500 uL
		20169-1mg	1 mg (lyophilized)
CF®543	541/560 nm	20316-1	50 uL
		20316	500 uL
		20316-1mg	1 mg (lyophilized)
CF®555	555/565 nm	20276-1	50 uL
		20276	500 uL
		20276-1mg	1 mg (lyophilized)
CF®568	562/583 nm	20377-50uL	50 uL
		20377-500uL	500 uL
		20377-1mg	1 mg (lyophilized)
CF®594	593/614 nm	20170-1	50 uL
		20170	500 uL
		20170-1mg	1 mg (lyophilized)
CF®633	630/650 nm	20171-1	50 uL
		20171	500 uL
		20171-1mg	1 mg (lyophilized)
CF®660C	667/685 nm	20372-50uL	50 uL
		20372-500uL	500 uL
		20372-1mg	1 mg (lyophilized)
CF®680	681/698 nm	20241-1	50 uL
		20241	250 uL
CF®770	770/797 nm	20242-1	50 uL
		20242	250 uL

Bv: Bovine; Ch: chicken; Gt: goat; Hs: horse; Hu: human; Ms: mouse;  
Rb: rabbit; Rt: rat; Sh: Sheep; SHm: Syrian hamster

### Concentration

#### Liquid format:

2 mg/mL in pH ~7.4 PBS containing 50% glycerol, 2 mg/ml bovine serum albumin (IgG-free and protease-free) and 0.05% sodium azide.

#### Lyophilized format (after reconstitution):

2 mg/mL in pH ~7.4 PBS containing 15 mg/mL bovine serum albumin (IgG-free and protease-free) and 20 mg/mL trehalose.

### Storage and Handling

Store at -20°C, protected from light. Product is stable for at least 6 months from date of receipt when stored as recommended. Liquid format antibodies contain 50% glycerol and will not freeze at -20°C.

Reconstitution (lyophilized format only): add 0.5 mL dH<sub>2</sub>O and mix gently to dissolve. Store at -20°C, protected from light. Aliquot to avoid freeze-thaw cycles. Alternatively, glycerol can be added to the antibody so that it will not freeze at -20°C: add 0.25 mL dH<sub>2</sub>O to the lyophilized antibody and mix gently to dissolve, then add 0.25 mL glycerol and mix well. Optional: a preservative may be added, such as 0.05% (final concentration) sodium azide.

**Note:** storage of the antibody for more than a day at final working dilution is not recommended.

### Product Description

Donkey anti-guinea pig IgG (H+L), whole antibody, highly cross-adsorbed is an affinity purified donkey antibody that recognizes guinea pig IgG. This antibody is highly-cross adsorbed for multiple labeling, and has minimal cross-reactivity with bovine, chicken, goat, Syrian hamster, horse, human, mouse, rabbit, rat and sheep serum proteins.

CF® Dyes are superior to Alexa Fluor®, Dylight®, and Cy® dyes for antibody labeling by having combined advantages in brightness, photostability, specificity and novel features ideal for in vivo imaging. For more information about CF® Dyes, download our CF® Dye Selection Guide at [www.biotium.com](http://www.biotium.com).

### Recommended Dilution Range

Fluorescence microscopy: 1-2 ug/mL

Flow cytometry: 1 ug/10<sup>6</sup> cells

Near-infrared western detection: 50-100 ng/mL

These concentrations are provided as a starting point for optimization; appropriate dilutions should be determined empirically. Generally IgG conjugates are used in the range of 1-10 ug/mL.

For more detailed protocols for immunofluorescence staining for fluorescence microscopy, flow cytometry, or fluorescence-based western detection, please visit our website.

## Accessory products for immunofluorescence

Cat. No.	Product	Features
23001	EverBrite™ Mounting Medium	<ul style="list-style-type: none"> <li>• Excellent protection from photobleaching</li> <li>• Compatible with a wide variety of dyes including cyanine dyes</li> <li>• Available in wet-set or hardset formulations</li> <li>• Available With or without DAPI</li> <li>• EverBrite TrueBlack® Hardset is the only mounting medium with autofluorescence quenching</li> <li>• Hardset and EverBrite TrueBlack® formulations available with NucSpot® 640 far-red nuclear counterstain</li> </ul>
23002	EverBrite™ Mounting Medium with DAPI	
23003	EverBrite™ Hardset Mounting Medium	
23004	EverBrite™ Hardset Mounting Medium with DAPI	
23015	EverBrite™ Mounting Medium with NucSpot® 640	
23016	EverBrite™ Hardset Mounting Medium with NucSpot® 640	
23017	EverBrite TrueBlack® Hardset Mounting Medium	
23018	EverBrite TrueBlack® Hardset Mounting Medium with DAPI	
23019	EverBrite TrueBlack® Hardset Mounting Medium with NucSpot® 640	
23005	CoverGrip™ Coverslip Sealant	<ul style="list-style-type: none"> <li>• Replaces nail polish for coverslip sealing</li> <li>• Won't mix with aqueous mounting media</li> </ul>
23007	TrueBlack® Lipofuscin Autofluorescence Quencher	<ul style="list-style-type: none"> <li>• Eliminates lipofuscin autofluorescence with less background than Sudan Black B</li> <li>• TrueBlack® Plus offers even lower far-red background than original TrueBlack®</li> <li>• Reduces background from other sources</li> <li>• Can be used before or after IF staining</li> </ul>
23014	TrueBlack® Plus Lipofuscin Autofluorescence Quencher	
40061	RedDot™2 Far Red Nuclear Counterstain	<ul style="list-style-type: none"> <li>• Far-red nuclear dye for the Cy®5 channel</li> <li>• More specific than Draq®7</li> </ul>
22005	Mini Super <sup>HT</sup> Pap Pen 2.5 mm tip, ~400 uses	<ul style="list-style-type: none"> <li>• Create hydrophobic barriers around tissue sections</li> <li>• Heat-stable to 120°C</li> <li>• Insoluble in aqueous buffers, detergents, alcohol and acetone; can be removed with xylene</li> </ul>
22006	Super <sup>HT</sup> Pap Pen 4 mm tip, ~800 uses	
22003	Mini-Cell Scrapers	<ul style="list-style-type: none"> <li>• Collect cells from 96-well to 24-well plates</li> <li>• Polypropylene, sterile, disposable</li> </ul>
23006	Flow Cytometry Fixation/Permeabilization Kit	<ul style="list-style-type: none"> <li>• Ready-to-use fix/perm buffers for intracellular staining</li> </ul>
22015	Fixation Buffer	<ul style="list-style-type: none"> <li>• Convenient buffers and blocking agents for immunofluorescence or western</li> </ul>
22016	Permeabilization Buffer	
22017	Permeabilization and Blocking Buffer	
22010	10X Fish Gelatin Blocking Agent	
22011	Fish Gelatin Powder	
22014	30% Bovine Serum Albumin Solution	
22002	Tween®-20	

Please visit [www.biotium.com](http://www.biotium.com) to view our full selection of products featuring bright and photostable fluorescent CF® Dyes, including Mix-n-Stain™ antibody labeling kits, primary antibody conjugates, streptavidin, phalloidin, and other bioconjugates, as well as conjugates of biotin, HRP, AP, R-PE, APC and PerCP.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.

CF dye technology is covered by US and international patents. Alexa Fluor is a registered trademark of Molecular Probes; CYDYE is a registered trademark of GE Healthcare; Draq7 is a registered trademark of Biostatus Ltd; DyLight is a registered trademark of Pierce Biotechnology; TWEEN is a registered trademark of Uniqema Americas LLC.