

## Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

## Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

www.biotium.com

Revised: January 22, 2019



# **Product Information**

### Rabbit Anti-Sheep IgG (H+L), Whole Antibody, Highly Cross-Adsorbed (min x Hu), CF® Dye Conjugates

| CF® Dye | Ex/Em      | Catalog No. | Size               |
|---------|------------|-------------|--------------------|
| CF®488A | 490/515 nm | 20172-1     | 50 uL              |
|         |            | 20172       | 500 uL             |
|         |            | 20172-1mg   | 1 mg (lyophilized) |
| CF®543  | 541/560 nm | 20323-1     | 50 uL              |
|         |            | 20323       | 500 uL             |
|         |            | 20323-1mg   | 1 mg (lyophilized) |
| CF®555  | 555/565 nm | 20173-1     | 50 uL              |
|         |            | 20173       | 500 uL             |
|         |            | 20173-1mg   | 1 mg (lyophilized) |
| CF®633  | 630/650 nm | 20174-1     | 50 uL              |
|         |            | 20174       | 500 uL             |
|         |            | 20174-1mg   | 1 mg (lyophilized) |

Hu: human

#### Concentration

Liquid format: 2 mg/mL in pH  $\sim$ 7.4 PBS containing 50% glycerol, 2 mg/ml bovine serum albumin (IgG-free and protease-free) and 0.05% sodium azide.

Lyophilized format (after reconstitution): 2 mg/mL in pH ~7.4 PBS containing 15 mg/mL bovine serum albumin (IgG-free and protease-free) and 20 mg/mL trehalose.

#### Storage and Handling

Store at -20°C, protected from light. Product is stable for at least 6 months from date of receipt when stored as recommended. Liquid format antibodies contain 50% glycerol and will not freeze at -20°C.

Reconstitution (lyophilized format only): add 0.5 mL dH<sub>2</sub>O and mix gently to dissolve. Store at -20°C, protected from light. Aliquot to avoid freeze-thaw cycles. Alternatively, glycerol can be added to the antibody so that it will not freeze at -20°C: add 0.25 mL dH<sub>2</sub>O to the lyophilized antibody and mix gently to dissolve, then add 0.25 mL glycerol and mix well. Optional: a preservative may be added, such as 0.05% (final concentration) sodium azide.

Note: storage of the antibody for more than a day at final working dilution is not recommended.

#### **Product Description**

Rabbit anti-sheep IgG (H+L), whole antibody, highly cross-adsorbed is an affinity purified rabbit antibody that recognizes sheep IgG. This antibody is highly-cross adsorbed for multiple labeling, and has minimal cross-reactivity with human serum proteins.

CF® dyes are superior to Alexa Fluor®, DyLight®, and Cy® dyes for antibody labeling by having combined advantages in brightness, photostability, specificity and novel features ideal for in vivo imaging. For more information about CF® dyes, download our CF® Dye Selection Guide at www.biotium.com.

#### **Recommended Dilution Range**

Fluorescence microscopy: 1-2 ug/mL Flow cytometry: 1 ug/10<sup>6</sup> cells Near-infrared western detection: 50-100 ng/mL

These concentrations are provided as a starting point for optimization; appropriate dilutions should be determined empirically. Generally IgG conjugates are used in the range of 1-10 ug/mL.

For more detailed protocols for immunofluorescence staining for fluorescence microscopy, flow cytometry, or fluorescence-based western detection, please visit our website.

#### **Related Products**

| Cat. No. | Product Name   |  |
|----------|--|--|
| 23001    | EverBrite™ Mounting Medium                             |  |
| 23002    | EverBrite™ Mounting Medium with DAPI                   |  |
| 23003    | EverBrite™ Hardset Mounting Medium                     |  |
| 23004    | EverBrite™ Hardset Mounting Medium with DAPI           |  |
| 23005    | CoverGrip™ Coverslip Sealant                           |  |
| 23007    | TrueBlack™ Lipofuscin Autofluorescence Quencher        |  |
| 40061    | RedDot™2 Far Red Nuclear Counterstain                  |  |
| 22005    | Mini Super <sup>н⊤</sup> Pap Pen 2.5 mm tip, ~400 uses |  |
| 22006    | Super <sup>н⊤</sup> Pap Pen 4 mm tip, ~800 uses        |  |
| 22003    | Mini-Cell Scrapers                                     |  |
| 23006    | Flow Cytometry Fixation/Permeabilization Kit           |  |
| 22015    | Fixation Buffer  |  |
| 22016    | Permeabilization Buffer                                |  |
| 22017    | Permeabilization and Blocking Buffer                   |  |
| 22010    | 10X Fish Gelatin Blocking Agent                        |  |
| 22011    | Fish Gelatin Powder                                    |  |
| 22014    | 30% Bovine Serum Albumin Solution                      |  |
| 22002    | Tween®-20  |  |

Please visit www.biotium.com to view our full selection of products featuring bright and photostable fluorescent CF® dyes, including Mix-n-Stain ™ antibody labeling kits, primary antibody conjugates, streptavidin, phalloidin, and other bioconjugates, as well as conjugates of biotin, HRP, AP, R-PE, APC and PerCP.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.

CF dye technology is covered by pending US and international patents. Alexa Fluor is a registered trademark of Molecular Probes; CYDYE is a registered trademark of GE Healthcare; Draq7 is a registered trademark of Biostatus Ltd; DyLight is a registered trademark of Pierce Biotechnology; TWEEN is a registered trademark of Uniqema Americas LLC.