

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



Revised: September 18, 2023

Product Information

130

250

140

260

TEV site (cleaved)

310 320
KGDTSGDYKK ALLLLCGEDD ENLYFQ

Annexin V

Product List

Product	Cat. No.	Unit Size
Annexin V (His Tag)	29088-100ug	100 ug
	29088	1 mg
Annexin V (No Tag)	29089-100ug	100 ug
	29089	1 mg
	29089-5mg	5 mg

Storage and Handling

Store at -20°C. Annexin V was lyophilized from solution containing 25 nM HEPES, 140 mM NaCl, 1 mM EDTA, pH 7.4 (0.2 um filtered). Product is stable for at least 6 months from date of receipt when stored as recommended.

Reconstitute lyophilized Annexin V with sterile deionized water for a stock solution of 1 mg/mL (buffer components listed above). Mix gently for 30 to 60 minutes at room temperature to dissolve, avoiding vigorous mixing or vortexing. Reconstituted Annexin V without BSA is stable for 1 week at 4°C, or up to 3 months at -20°C under sterile conditions. Aliquots should not be refrozen after thawing.

For longer storage, BSA can be added to a final concentration of 1 mg/mL to stabilize Annexin V during freezing. Annexin V stock solutions with 1 mg/mL BSA can be stored for at least 6 months at 4°C as sterile solutions, or with the addition of 0.1% sodium azide as a preservative.

Note: Do not add BSA if conjugating this protein to a dye or other molecule.

Molecular Weight: 37.6 kDa (His Tag), 36.7 kDa (No Tag) **Purity:** >98% pure as demonstrated by reducing SDS-PAGE

Endotoxin Level: < 1 EU/ug, determined by PyroGene® Recombinant

Factor C Assay.

Product Description

Annexin V is a member of the annexin family of calcium-dependent phospholipid-binding proteins with high affinity for phosphatidylserine (PS). In viable cells, PS is predominantly located on the inner leaflet of the plasma membrane. During apoptosis, PS translocates from the inner leaflet to the outer leaflet. Annexin V can bind to PS located on the outer leaflet and therefore fluorescently-labeled Annexin V is commonly used as a readout of apoptosis, as well as a probe for staining extracellular vesicles (EVs). Annexin V has also been identified as an anticoagulant protein that is involved in the blood coagulation cascade through indirect inhibition of the thromboplastin-specific complex.

Recombinant Human Annexin V is expressed from *E. coli* as a single, non-glycosylated polypeptide chain. Our Annexin V is offered with and without a C-terminal 6X His tag.

Annexin V (His Tag) 20 30 40 MAQVLRGTVT DFPGFDERAD AETLRKAMKG LGTDEESILT LLTSRSNAQR QEISAAFKTL 100 FGRDLLDDLK SELTGKFEKL IVALMKPSRL YDAYELKHAL KGAGTNEKVL TEIIASRTPE 150 ELRAIKQVYE EEYGSSLEDD VVGDTSGYYQ RMLVVLLQAN RDPDAGIDEA QVEQDAQALF 190 200 210 220 230 QAGELKWGTD EEKFITIFGT RSVSHLRKVF DKYMTISGFQ IEETIDRETS GNLEQLLLAV 260 270 280 290 VKSIRSIPAY LAETLYYAMK GAGTDDHTLI RVMVSRSEID LFNIRKEFRK NFATSLYSMI 310 320 330 KGDTSGDYKK ALLLLCGEDD ENLYFQSHHH HHH TEV site Annexin V (No Tag) 30 MAQVLRGTVT DFPGFDERAD AETLRKAMKG LGTDEESILT LLTSRSNAQR QEISAAFKTL FGRDLLDDLK SELTGKFEKL IVALMKPSRL YDAYELKHAL KGAGTNEKVL TEIIASRTPE

Figure 1. Amino acid sequences for recombinant Annexin V (His Tag), top; and Annexin V (No Tag), bottom. C-terminal TEV protease cleavage site (blue) and 6X His tag (red) are indicated.

150

ELRAIKQVYE EEYGSSLEDD VVGDTSGYYQ RMLVVLLQAN RDPDAGIDEA QVEQDA

270

OAGELKWGTD EEKFITIFGT RSVSHLRKVF DKYMTISGFO IEETIDRETS GNLEOLLLAV

VKSIRSIPAY LAETLYYAMK GAGTDDHTLI RVMVSRSEID LFNIRKEFRK NFATSLYSMI

160

280

170

290

180

Related Products

Tolatou i Toudoto		
Cat. No	Product	
2900329085	Annexin V CF® Dye Conjugates	
2900629082	Annexin V Near IR CF® Dye Conjugates	
29003R 29069R	Annexin V CF® Dye Conjugates, Azide-Free, Lyophilized	
99902	5X Annexin Binding Buffer (15 mL)	
30119-30122	ExoBrite™ Annexin EV Staining Kits	
1040210408	NucView® Caspase-3 Enzyme Substrates	
30067	Dual Apoptosis Assay Kit with NucView® 488 Caspase-3 Substrate & CF®594 Annexin V	
30076	Dual Apoptosis Assay Kit with NucView® 488 Caspase-3 Substrate & CF®640R Annexin V	
30062	NucView® 488 and MitoView™ 633 Apoptosis Kit	
30072	NucView® 488 and RedDot™2 Apoptosis and Necrosis Kit	
30065	Apoptosis & Necrosis Quantitation Kit Plus	
30066	Apoptotic, Necrotic & Healthy Cells Quantitation Kit Plus	
30060	CF®488A Annexin V and 7-AAD Apoptosis Kit	
30061	CF®488A Annexin V and PI Apoptosis Kit	
40061	RedDot™2 Far-Red Nuclear Stain (for dead or fixed cells)	
4008341040	NucSpot® Nuclear Stains (for dead or fixed cells)	
40085	NucSpot® Far-Red, 1000X in DMSO (for dead or fixed cells)	
40084	7-AAD Solution, 1 mg/mL	
30001	JC-1 Mitochondrial Membrane Detection Kit	
70076	Aquaphile™ JC-1	
30063	CF®488A TUNEL Assay Apoptosis Detection Kit	
30064	CF®594 TUNEL Assay Apoptosis Detection Kit	
30074	CF®640R TUNEL Assay Apoptosis Detection Kit	
32010	Live-or-Dye NucFix™ Red Staining Kit	
32002-32018	Live-or-Dye™ Fixable Viability Staining Kits	
00025	Staurosporine	

Please visit our website at www.biotium.com for information on our life science research products, including probes for monitoring cell viability, cellular stains, Annexin V conjugates, and other bioconjugates.

Pyrogene is a registered trademark of Lonza Bioscience.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.